



PRICON

we share healthcare

PRODUCT CATALOUGE

Ophthalmics



MICRO SURGICAL OPHTHALMIC CANNULAE & INSTRUMENTS

Having aggressively addressed ourselves to the myriads of the fast developing microsurgical needs,  today is widely referred as THE SOURCE for Ophthalmic Cannulae and instruments having exports to many countries including USA, UK, Europe, Far East, Middle East, Africa and Russia.

CONSISTENT QUALITY

We are proud to admit that Quality Control and specification adherence are the hallmark of the Company. As a leader in Quality and Innovation, We are committed to distancing itself from our competitors by supplying the high quality products at competitive prices without compromising in materials and workmanship. For the Customer, our benefits are obvious. Quite simply  Products perform better because they are made better.



ISO 9001:2015



EN ISO 13485:2016



CE



WHO-GMP

THE R & D

Products are just example of our commitment to the profession of surgery. Our R&D Centre is the core of an ongoing research programme in consultation with the global surgeons to produce more efficient and more variety of instruments. We welcome and would like to hear from you about your thoughts and ideas on new designs on improvements to our current line.

We welcome your new product idea which can simplify surgical techniques and reduce procedural costs.

This is our representative catalogue, for further details please visit our website.

REUSABLE & DISPOSABLE CANNULAE

ANESTHESIA NEEDLES
CYSTOTOMES
LACRIMAL
AIR INJECTION & IRRIGATING CANNULAE
HYDRODISSECTION / HYDRODELINEATION
CORTEX ASPIRATION
IRRIGATING & ASPIRATING
CAPSULE POLISHERS
LENS REMOVAL CANNULAE
VITREORETINAL
BIMANUAL ASPIRATION / IRRIGATION
LASIK
GLAUCOMA
IRRIGATING CHOPPER
MISCELLANEOUS & NEW ARRIVALS

MICROSURGICAL INSTRUMENTS (STEEL / TITANIUM)

SPECULUMS
RETRACTORS
LACRIMAL & DCR INSTRUMENTS
HOOKS
SPATULAS
LOOPS
SPOONS & CURETTES
NUCLEUS ROTATORS / MANIPULATORS
FIXATION RINGS
CALIPERS / GAUGES
FORCEPS
SCISSORS
NEEDLE HOLDERS
BLADE BREAKERS / HOLDERS
STERILIZATION TRAYS
MARKERS
DALK
DSEK

SINGLE USE INSTRUMENTS

MANIPULATORS

CHOPPER

MISCELLANEOUS INSTRUMENTS

IRIS RETRACTOR

EYE SHIELD

SPONGE

TREPHINES

DIAGNOSTIC STRIPS

ARTIFICIAL IMPLANTS

TROCAR SET

DISPOSABLE KNIVES

PRICON GOLD - JGH

SINGLE USE WITH HANDLE - JSH

SINGLE USE WITH HUB - JS

SURGERY SETS RANGE OF PRODUCTS

INSTRUCTION FOR ORDERING

CANNULAE:

The Cannulae are manufactured both in reusable & disposable (Single use). While ordering for single use add "D" after the catalogue number.

INSTRUMENTS

The instruments are manufactured in stainless steel & titanium .
While ordering please add "T" for titanium instruments & "D" for disposables.

BLADES

The microsurgical blades are manufactured in three types :-

- | | |
|---|--------|
| (a) Pricon Gold with handle (sterile) & Chemically Sharpened Blades | -- JGH |
| (b) Single use with handle (sterile) Standard Type | -- JSH |
| (c) Single use with hubs (sterile) Economy Type | -- JS |

While ordering add all suitable abbreviations after catalogue number referred above.


PACKING : The Cannulae & Blades are packed in box. The instruments are packed individually.

STERILITY : All single use products are ETO sterile.

Anaesthesia Needles

RETROBULBAR NEEDLE (ATKINSON)

[R] [S]



JD001 - 27	27G x 38mm	JD001.1 - 27	27G x 32mm
JD001 - 25	25G x 38mm	JD001.1 - 25	25G x 32mm
JD001 - 23	23G x 38mm	JD001.1 - 23	23G x 32mm

BOX OF 10

RETROBULBAR NEEDLE (SHARP POINT)


[R] [S]



JD003 - 20	20G x 38mm
JD003 - 23	23G x 38mm
JD003 - 25	25G x 38mm

PERIBULBAR NEEDLE (ATKINSON)

[R] [S]



JD005 - 27	27G x 22mm	JD005.1 - 27	27G x 32mm
JD005 - 25	25G x 22mm	JD005.1 - 25	25G x 32mm
JD005 - 23	23G x 22mm	JD005.1 - 23	23G x 32mm

SUBTENON'S ANESTHESIA

[R] [S]



JD009 - 19	19G x 25mm
------------	------------

Curved Flattened Tip

SUBTENON'S ANESTHESIA

[R] [S]



JD011C - 21	21G x 25mm	JD011A - 21	21G x 25mm
JD011C - 19	19G x 25mm	JD011A - 19	19G x 25mm
Curved		Angled	

SUBTENON'S ANESTHESIA PARA NEEDLE

[R] [S]



JD009-19

Curved

SUBTENON'S ANESTHESIA - CURVED

[R]



J1009 - 23	23G x 35mm
------------	------------

Curved Flattened Bulbous Tip, Top Port

MUTHUSAMY ANESTHESIA - NEEDLE

[R] [S]




JD1013 - 24	24G x 32mm
-------------	------------

Design facilitates to course around the globe and parallel to the optic nerve.

Cystotomes

CYSTOTOME (STRAIGHT)


[S]



Jd015 - 25	25G x 16mm
JD015 - 27	27G x 16mm
JD015 - 30	30G x 16mm

IRRIGATING CYSTOTOME (FORMED)


[S]



JD017 - 25	25G x 12mm	JD017.1 - 25	25G x 8.5mm
JD017 - 26	26G x 12mm	JD017.1 - 26	26G x 8.5mm
JD017 - 27	27G x 12mm	JD017.1 - 27	27G x 8.5mm
JD017 - 30	30G x 12mm	JD017.1 - 30	30G x 8.5mm

CYSTOTOME (ANGLED)

[S]




JD019 - 25	25G
JD019 - 27	27G
JD019 - 30	30G

Tip Rotated 45°

IRRIGATING CYSTOTOME (REVERSE)

[S]



JD021 - 21	21G x 22mm
------------	------------

Tip Bent towards Bevel

SHARP IRRIGATING CYSTOTOME (KELMAN)

[S]




JD023 - 25	25G x 32mm
------------	------------

Curved Sharp Cutting Edge

IRRIGATING CYSTOTOME (BLUMENTHAL)

[S]



JD025 - 25	25G x 22mm
------------	------------


Formed 8.5mm from bent to tip

Lacrimal

LACRIMAL CANNULA REINFORCED (STRAIGHT) R S



JD027- 23 23G x 32mm
JD027- 25 25G x 32mm



20G Shaft having reinforced
12mm 23G cannula (Malleable), Straight

LACRIMAL CANNULA REINFORCED (CURVED) R S



JD029-23 23G x 32mm
JD029-25 25G x 32mm




20G Shaft having reinforced
12mm 23G cannula (Malleable), Curved

LACRIMAL CANNULA (BAILEY) R S



JD031- 23 23G x 15mm
JD031- 25 25G x 15mm




For Infants short 20G Straight Shaft
having 7mm Exposed Tip of 23G

LACRIMAL CANNULA (SHAHINIAN) R S




JD1031- 23 23G x 13mm JD1031.1 - 23 23G x 35mm
JD1031- 25 25G x 13mm JD1031.1 - 25 25G x 35mm



LACRIMAL CANNULA (FASANELLA) R S



JD033-23 23G x 28mm
JD033-25 25G x 28mm
JD033-26 26G x 28mm



JD033.1-23 23G x 42mm
JD033.1-26 25G x 42mm

LACRIMAL CANNULA (ANEL) BLUNT TIP R S



JD035S - 23
JD035C - 23



23G Straight
23G Curved



JD035S - 25 25G Straight
JD035C - 25 25G Curved

LACRIMAL CANNULA (WEST) R S



JD037-23 23G x 30mm



Gently Curved 0.3mm Bottom Port


FRONTALIS SUSPENSION SET




JD041- 20G



Silicone Sleeve



Stainless steel needles 0.90 mm x 6.3 cm Solid silicone tubing length 40 cm diameter 0.80 mm.
Silicone tubing sleeve length 10 mm ID - 1 mm; OD - 2 mm



LACRIMAL INTUBATION SYSTEM



LACRIMAL INTUBATION SET



JD043 - 23	23G/11cm Probe	JD043.2-23	23G/17.5cm Probe
JD043 - 25	25G/11cm Probe	JD043.2-25	25G/17.5cm Probe
JD043 - 27	27G/11cm Probe	JD043.2-27	27G/17.5cm Probe
Stainless Steel probes with 30cms.			
JD043.1-23	23G/7.5cm Probe		
JD043.1-25	25G/7.5cm Probe		
JD043.1-27	27G/7.5cm Probe		


DACRYOCYSTORHINOSTOMY (DCR)



JD045S-20G	Length 4.50cms.	JD045A-20G	Length 4.50cms.
JD045S-23G	Length 4.50cms.	JD045A-23G	Length 4.50cms.
JD045.1S-20G	Length 3.5cms.	JD045.1A-20G	Length 3.5cms.
JD045.1S-23G	Length 3.5cms.	JD045.1A-23G	Length 3.5cms.

Stainless Steel Probes solid Silicone Rod


LACRIMAL INTUBATION SET (OLIVE TIP)



JD1041-23	23G	JD1041-23RD	23G
JD1041-25	25G	JD1041-25RD	25G
JD1041-27	27G	JD1041-27RD	27G

Stainless Steel probes 0.40mm x 11cms with Olive Tips. Silicone Tube length 35cms.

LACRIMAL INTUBATION SET (NOTCHED)



JD1043-23	23G
JD1043RD-23	23G

Stainless Steel probes 0.60mm x 11cms with Notched Tips. Silicone Tube Length 35cms.


LACRIMAL INTUBATION SET WITH SUTURE (CRAWFORD)



JD1045-20 23G

Stainless Steel probes 0.60mm x 11cms with Olive Tips, 6.0 Silk thread inserted inside the lumen of probe for easy removal without anesthesia


LACRIMAL INTUBATION SET (WITH SHEATH)




JD047-21	21G	JD048-21	21G	(Adult) Stainless Steel probes 0.50mm x 21cms
JD047-23	23G	JD048-23	23G	with Sheath Silicone Tube Length 40cms.
JD047-25	25G	JD048-25	25G	(Child) Stainless Steel probes 0.50mm x 18cms
				with Sheath Silicone Tube Length 40cms.

ANTERIOR CHAMBER/VISCOELASTIC IRRIGATING CANNULA


Anterior Chamber Cannula Rycroft

		4MM Bend	5MM Bend
	21G	JD051-21	JD051.1-21
	22G	JD051-22	JD051.1-22
	23G	JD051-23	JD051.1-23
	25G	JD051-25	JD051.1-25
	26G	JD051-26	JD051.1-26
	27G	JD051-27	JD051.1-27
	30G	JD051-30	JD051.1-30

Anterior Chamber Cannula Angular

		6MM Bend	7MM Bend	8MM Bend
	21G	JD052-21	JD052.1-21	JD052.2-21
	22G	JD052-22	JD052.1-22	JD052.2-22
	23G	JD052-23	JD052.1-23	JD052.2-23
	25G	JD052-25	JD052.1-25	JD052.2-25
	26G	JD052-26	JD052.1-26	JD052.2-26
	27G	JD052-27	JD052.1-27	JD052.2-27
	30G	JD052-30	JD052.1-30	JD052.2-30

Viscoelastic Irrigating Cannula

		9MM Bend	10MM Bend	11MM Bend	12MM Bend
	21G	JD053-21	JD053.1-21	JD053.2-21	JD053.3-21
	22G	JD053-22	JD053.1-22	JD053.2-22	JD053.3-22
	23G	JD053-23	JD053.1-23	JD053.2-23	JD053.3-23
	25G	JD053-25	JD053.1-25	JD053.2-25	JD053.3-25
	26G	JD053-26	JD053.1-26	JD053.2-26	JD053.3-26
	27G	JD053-27	JD053.1-27	JD053.2-27	JD053.3-27
	30G	JD053-30	JD053.1-30	JD053.2-30	JD053.3-30

Reusable & Disposable Cannulae



Air Injection & Irrigating Cannulae

IRRIGATING CANNULA STRAIGHT


R S



JD049-20 20Gx25mm JD049-23 23Gx25mm
JD049-22 22Gx25mm JD049-25 25Gx25mm
Blunt Tip

VISCO FLOW CANNULA

R S



JD1056-27 27G x 24mm

For easier injection of viscoelastic to all meridians Dome shape provides safe insertion.

BISHOP HARMON CANNULA

R S




JD055-16 16G x 25mm
JD055-19 19G x 25mm

40° Angled, Spatulated, 8mm from end

ANTERIOR CHAMBER CANNULA (MCINTYRE)

R S



JD061 - Set	Set of 5 Cannulae	JD061 - 26	26G Angled
JD061 - 30	30G Angled	JD061 - 25	25G Angled
JD061 - 27	27G Angled	JD061 - 23	23G Angled

FLATTENED ANTERIOR CHAMBER CANNULA (BRACKEN)

R S



JD057-19 19G x 23mm

Curved 90° flattened bevel tip

ANTERIOR CHAMBER MAINTAINER (LEWICKY)

R S



JD065-20 20G

Silicone tubing attached to 20G x 3.50mm cannulae unique annular rings for self retaining.

ANTERIOR CHAMBER CANNULA (KNOLLE)

R S



JD063 - 23	23G x 19mm	JD063-Set	Set of 3 Cannulae
JD063 - 25	25G x 19mm		
JD063 - 30	30G x 19mm		2mm Angled Tip

ANTERIOR CHAMBER MAINTAINER BULBOUS TIP (JAYMCE)

R S




JD1067-20 20G x 4mm

Bulbous tip allows easy insertion and minimum tissue disturbance

ANTERIOR CHAMBER MAINTAINER (BLUEMENTHAL)

R S



JD1065-23 23G x 05mm
20cms Silicone tubing attached to reinforced Cannula with Atkinson point. Slightly spatulated for self retaining.

NUCLEUS HYDRODISSECTOR (RAININ) HORIZONTAL

R S



JD069H-23 23Gx22mm	JD069H-27 27Gx22mm
JD069H-25 25Gx22mm	JD069H-30 30Gx22mm

Angled 45°, Bend 7mm Spatulated on horizontal Plane.

Hydrodissection / Hydrodelineation

NUCLEUS HYDRODISSECTOR (RAININ) VERTICAL


R S



JD069V-27 27Gx22mm
JD069V-25 25Gx22mm
Angled 45°, Bend 7mm Spatulated on Vertical Plane.

NUCLEUS HYDRODISSECTOR (MCINTYRE)

R S



JD071-23 23Gx22mm
JD071-25 25Gx22mm
JD071-27 27Gx22mm
JD071-30 30Gx22mm

Angled 45°, Bend 11mm Flattened on horizontal Plane.

HYDRODISSECTOR (PEARCE)

R S



JD1073 - 25G
JD1073 - 27G
Micro size spatulated tip for easy entry through wound under capsule Hook 2.00mm x 2.50mm

MICRO HYDRODISSECTOR (PEARCE)

R S



JD073-25 25Gx22mm
JD073-27 27Gx22mm
JD073-30 30Gx22mm


Hook 1.75mm wide 2.25mm long length 22mm

Reusable & Disposable Cannulae



Hydrodissection / Hydrodelineation

NUCLEOUS HYDROLYSIS NEEDLE (P.KOCH) R S



JD075-25 25G x 22 mm




Curved Short bevel tip penetrates the nucleus for manual sculpting to improve lamellar separation. Flattened & Smooth insertion under anterior capsule.

NUCLEOUS HYDRODISSECTOR (PEARCE) R S



JD1075-25 25G x 22mm




Close tip to prevent clogging while entry. Three Irrigating ports for improved fluid distribution and lamellar separation.

NUCLEOUS HYDRODelineATOR (FEASTER) R S



JD1077-25 25G x 28mm




Closed wedge shaped tip allows easy penetration between nucleus & epinucleus. Two irrigating ports for fluid distribution in both the planes

NUCLEOUS HYDRODelineATOR (BLUMENTHAL) R S



JD1078-25 25G x 22mm




Tapered tip allows more flow of fluid through a smaller tip. Easy insertion under capsule for hydrodissection or into nucleus for hydrodelineation.

NUCLEOUS HYDRODISSECTOR R S



JD1079-27 27G x 22mm




Horizontally Flattened Tip

HYDRODISSECTOR CANNULA (CHANG) R S



JD1080-27 27G




27G, 90° angle, flat tip with a beveled opening. Flat tip for easy insertion under the capsular rim, provide broad stream of fluid Beveled tip used to rotate the nucleus to ensure that it is completely separated.

Cortex Aspiration

CORTEX ASPIRATING CANNULAE (MCINTYRE BINKHORST) R S



JD079-R-26 26G x 22mm
JD079-L-26 26G x 22mm




Blunt Tip, Front opening, Hooked Left and Right

CORTEX ASPIRATOR (ALPERN) R S



JD081-R-25 25G x 22mm
JD081-L-25 25G x 22mm




Angled 450,5 mm from end
Hook 1.50 mmx1.75 mm

ASPIRATING CANNULA (GILLS) R S



JD085-5 5mm JD085-8 8mm
JD085-6 6mm JD085-10 10mm




Angled to tip

ASPIRATION CANNULA (HEALON) R S



JD087- 19 19Gx32mm
JD087- 21 21Gx32mm
JD087- 22 22Gx32mm




Angled 10 mm from tip

CORTEX EXTRACTOR CANNULA (SIMCOE) R S



JD091-21 21Gx21mm JD091.1-21 21Gx13mm
JD091-23 23Gx21mm JD091.1-23 23Gx13mm




Smooth closed tip with top aspirating port into curve

CORTEX EXTRACTOR CANNULA R S



JD092-19 19G




W/Adapter, W/O Polisher Simcoe

CORTEX EXTRACTOR CANNULA (WELSH) R




J095-25 25G
J095-S 25G Straight J095-R 25G Right
J095-L 25 G Left Flat Olive Tip length 25 mm



ALPHA CHROMOSTRYPsin ASPIRATING CANNULA (TROUTMAN) R



J097-25 25G x 26mm



Smooth bulbous tip

Reusable & Disposable Cannulae



Cortex Aspiration

U SHAPED ASPIRATING CANNULA

R S





JD099-21 21Gx22mm JD099-25 25Gx22mm
JD099-23 23Gx22mm JD099-27 27Gx22mm

Hooked 1.50 mmx1.75 mm wide

CORTEX EXTRACTOR (CHARLEUX)

R S




JD1098-24 24G

Bulbous tip allows easy insertion and minimum tissue disturbance.

Irrigating & Aspiring

DOUBLE I & CANNULA


R



J103-21 21G
J103-21 21G Thinwall

I & CANNULA (GILLS)

R




J105-23 23G
J105-25 25G Thinwall

Angled 35, 8 mm from bend to tip

J SHAPED I & CANNULA

R




J107-23 23G

23G/23G Open end Aspiring Port

J SHAPED I & CANNULA (INSIDE PORT)

R




J109-23 23G

23G/23G Aspiring cannula has 0.30 mm inside port

CO-AXIAL I & A CANNULA (PEARCE)

R




J113-C Curved
J113-A Angled
J113-S Straight

18 G Outer wall 23 G Inner with 0.30 mm side Aspiring Port

BASIC COAXIAL CANNULA (MCINTYRE)

R




J115 0.3mm Port
J115-A 0.2mm Port

With Nylon connector & Silicone tubing

I & A CANNULA (GILLS)

R




J119-23 23G
J119-23 23G Reverse

Side by side front opening cannula with silicone tubing

I & A CANNULA (GILLS/WELSH)

R




J121-25 25G
J121-S Straight J121-L Left J121-R Right
J123 Same but Reverse Model

Flat Olive Tip, 25 G, Aspiration through hub 23 G Irrigation with Silicone tubing

I & CANNULA (SIMCOE REGULAR)

R S




JD129-20 20G 0.50mm Port
JD129-21 21G 0.45mm Port
JD129-22 22G 0.40mm Port
JD129-23 23G 0.30mm Port
JD129-25 25G 0.25mm Port

Aspiration on Right for Left handers

I & A CANNULA (SIMCOE REVERSE)

R S



JD133-22 22G 0.40mm Port
JD133-23 23G 0.30mm Port
JD133-25 25G 0.25mm Port

Aspiration through hub 0.30 mm Port, Irrigation through tubing

Reusable & Disposable Cannulae



Irrigating / Aspirating

I & A CANNULA (SIMCOE-DAHAN)


R S



JD135
JD137 Reverse Model

JD135-22 22G 0.40mm Port JD137-22 22G 0.40mm Port
JD135-23 23G 0.40mm Port JD137-22 23G 0.40mm Port (Reverse Model)

Same as JD129 but both the cannula tips are sand blasted at bottom to facilitate polishing poster or anterior capsules




"J" SHAPED I & A CANNULA (REGULAR)

R S




JD1149L-23 23G
JD1149R-23 23G

"J" shape facilitates removal of cortex at 120' clock position
JD1149 Regular JD1149-L Left JD1149-R Right
JD1150 Reverse JD1150-L Left JD1150-R Right 0.30 mm top port



WB I & A CANNULA (SIMCOE) WITH SILICONE BULB


R



J131


J131 WB-3 0.30 mm poprt J131 WB-4 0.40 mm port

Specially designed hub eliminates use of adpoter and directly fits in the silicon bulb. Inside taper 6% luer, can also be used with syringe. Regular.



I & A CANNULA (MODIFIED) (SIMCOE)


R



J1151-23 23G
J1152-23 23G

J1152 Same as 1151 but Reverse Model.

Aspiration through 23G tubing and irrigation through 0.50 mm (22 G) bottom port, Bottom irrigation keeps capsules away from the aspiration port.



Capsule Polishers


CAPSULE POLISHER (ANIS)

R S



JD155-S-23 23G x 22mm
JD155-A-23 23G x 22mm

Sand blasted 0.3 mm top irrigating port




CAPSULE POLISHER (JENSON)

R S



JD157-25 25G
JD157-27 27G

All around sand blasted bulbous tip, Length 28 mm, gently curved Bottom sand blasted




OLIVE TIP POLISHER (PEARCE) (TOP PORT)

R S



JD1157-25 25G x 25mm

Olive shaped tips with 0.30 mm top port micro sand blasted




POSTERIOR CAPSULE POLISHER (FREEMAN)

R S



JD159-19 19G x 22mm

Irrigating port on both sides Angled, Sand blasted




POSTERIOR CAPSULE POLISHER (SIMCOE)

R S



JD161-23 23G x 16mm

Reverse 0.30 mm aperture, Curved, Sand blasted




POSTERIOR CAPSULE POLISHER (SIMCOE)

R S




JD163-19 19Gx16mm
JD163-21 21Gx16mm
JD163-23 23Gx16mm

Top port 0.30 mm aperture curved, Sand blasted




ANTERIOR CAPSULE SCRUBBER (SIMCOE)

R S




JD165-23 23G x 19mm

Tip is Angled 45°, and roughened slightly




POSTERIOR CAPSULE SCRUBBER (SIMCOE)

R S



JD167

Same as JD165 but for posterior capsule



Reusable & Disposable Cannulae


Capsule Polishers

CAPSULE POLISHER (KRATZ)

RS




JD169-21 21Gx22mm
JD169-23 23Gx22mm



0.30 mm top port, Angled, Blasted 3 mm

COLLAR BUTTON (GREATHER)

RS



JD171-S-23 23G x 25mm
JD171-A-23 23G x 25mm




Micro Iris retractor sand blasted posterior surface for capsule polishing side irrigation

FLEXIBLE TIP POLISHER (ENGLELS)

RS



JD173-23 23Gx22mm
JD173-25 25Gx22mm
JD173-27 27Gx22mm




Curved, 23mm sleeve

FLEXIBLE TIP POLISHER

RS



JD175-23 23Gx12mm
JD175-25 25Gx12mm
JD175-27 27Gx12mm



Silicon Sleeve Anterior Capsule Polisher
Angled 45, bend 9mm

Lens Removal Cannulae

IRRIGATING VECTUS (KNOLE-PEARCE)

RS



JD183-23 23G Three Port
JD183-25 25G Two Port
JD183-27 27G Single Port



Loop 5 mm Wide, 9 mm Long

REVERSE IRRIGATING VECTUS (DREWS-KNOLLE)

RS



JD185-23 23G



One forward & two reverse irrigating ports,
Loop 6 mm, 9 mm Long, Overall length 35 mm

IRRIGATING VECTUS (SHEETS)

RS



JD187-21 21G x 19mm




3 Irrigating ports, Serrations on top of
vectus 21 G Loop 6.5 mm Wide,
19 mm Long, Overall length 38 mm

IRRIGATING VECTUS (ANIS)

RS



JD189-23 23G




Two front opening ports, 23 Gauge heart
shaped loop

IRRIGATING VECTUS (DREWS)

RS



JD191-25 25G



Posterior capsule polisher lightly angled
Sand blasted base irrigating port
facing inside circle 25 G.


NUCLEUS DELIVERY CANNULA (SIMCOE)

R



J193

J193-L Left J193-R Right with serrations




IRRIGATING VECTUS (MICRO)

RS



JD195-25 25G



Lightly Angled, Bottom sand blasted,
irrigating port facing outside circle.

SMALL INCISION NUCLEUS DELIVERY CANNULA WITH SPOON (MCINTYRE)

R



JD197-23 23G x 2.5mm Spoon




Used with trisector to take out bigger fragments.
3 mm / 4.5 mm width 23 G.

NUCLEUS DELIVERY PUNCH (BIDAYE)

RS



JD199-23 23G




The ridge on the loop entangles nucleus
when punched into it and delivers
it out of ECC inot A/C. Tapered plate prevents
Iris prolapse. J199B With Mesh

FISH HOOK CANNULA

RS



JD200-30 30G



J shaped Sharp Tip for sophisticated nucleus

Reusable & Disposable Cannulae


Vitreoretinal

HIGH VISCOSITY INJECTOR (PEYMAN)

R S



JD251-16 16Gx09mm
JD251-19 19Gx09mm



To facilitate inject aspirate of silicone oil in the repair of retinal detachment

HIGH VISCOSITY INJECTOR 4 mm (HAMMER)

R S



JD253-20 20G 2.5mm
JD253.1-20 20G 04mm
JD253.2-20 20G 06mm



Tubing 16 cms long designed for use with silicone oil syringe. 20G ultra thin wall cannula with 45° bevel tip

INFUSION CANNULA

R S



JD257-20 20G 2.5mm JD257.1-20 20G 4.0mm JD257.2-20 20G 6.0mm
JD257-23 23G 2.5mm JD257.1-23 23G 4.0mm JD257.2-23 23G 6.0mm
JD257-25 25G 2.5mm JD257.1-25 25G 4.0mm JD257.2-25 25G 6.0mm




SUBRETINAL INJECTOR (GLASER)

R S



JD263
JD263-S Straight
JD263-C Curved




20G x 25 mm Cannula with 32G (0.32 mm) x 3 mm tip extension

SUBRETINAL ASPIRATOR (LANDERS)

R S



JD267




20G x 25 mm Cannula with 27G x 6 mm tip extension Angled 60°, tip to aspirate lysed blood cells and subretinal fluid.

VITREORETINAL CANNULA (VANDER)

R S



JD269
JD269-S Straight
JD269-A Angled 45°



20G x 1" Cannula with 30 G X 5 mm tip

VITREORETINAL ASPIRATOR (MAUMENEE)

R S



JD273




20G Ultra thin wall straight shaft with graduation at 5 mm interval

SUBRETINAL FLUID CANNULA 6 mm (GRIZZARD)

R S



JD275
JD275-S Straight
JD275-C Curved




20Gx22 mm silicone extends 6 mm beyond end. Overall length 35 mm flexible tip allows automatic entry through retinal tears.

SILICON TIP CANNULA

R S



JD277-23 23G
JD277-25 25G




SUBRETINAL FLUID CANNULA 3 mm (GRIZZARD)

R S



JD279S-20 20G
JD279C-20 20G



20G x 22 mm Sleeve extends 3 mm beyond end overall length with tubing

FLUTE CANNULA (CHARLES) WITH VACUUM HANDLE

R



J283-20 20G x 25mm



Vacuum handle with 20 G x 25 mm cannula provides an exit port for blood or debris when used with independent infusion

VITREORETINAL BRUSH (GRIZZARD)

R S



JD285C-20 20G
JD285S-20 20G
JD285S-20 20G



Sleeve extends 3 mm beyond end have slits for atraumatic brushing of retina

PERFLUOROCARBON CO-AXIAL I/A CANNULA

R S



JD289-19 19G
JD289-23 23G
JD289-25 25G




Used for perfluorocarbon liquid exchange during retinal procedures.

PFCL CANNULA (NAGPAL)

R



J291



26G/26G Thinwall, Angled 35, 8 mm from bend to tip

Reusable & Disposable Cannulae

Vitreoretinal

MICROPICK GLASER 25G, 32 MM



JD293-23
JD293-25
JD293-27



PFCL EXCHANGE MULTI BORE CANNULA, 23G




JD511-23
JD511-25



VFI CANNULA (POLY TIP)



JD515-23
JD515-25
JD515-27



SILICON OIL REMOVAL CANNULA




JD519-23
JD519-25
JD519-27



Flut Cannula




JD500-20 20G JD500-23 23G
JD500-22 22G JD500-25 25G





SUPRACHOROIDAL INJECTION NEEDLE 30G

The suprachoroidal space (SCS) is a potential space between the sclera and choroid that traverses the circumference of the posterior segment of the eye. The SCS is an attractive site for drug delivery because it targets the choroid, retinal pigment epithelium, and retina with high bioavailability, while maintaining low levels elsewhere in the eye.


1000 Microns Out
of the 30G needle



JD530-30G



Before Injection



SCS Injection

Suprachoroidal
space

Product Code : JD530-30G Suprachoroidal Injection Needle

Reusable & Disposable Cannulae


Lasik

DALK AIR INJECTION

[R]



JD361-27 27G x 32mm



Big Bubble 27 G Air Injection Cannula


LASIK IRRIGATING CANNULA (BANAJI)

[R] [S]



JD373-25 25G
JD373-27 27G

Design permits sideways entry into the flap without dragging epithelium the interface. Irrigation by four 0.25 mm ports on periphery




LASIK FLAP IRRIGATOR (VIDAURI)

[R] [S]



JD371-16 16G x 12mm

Universal cannula simultaneously irrigates both sides of flap. Washes at Stromal bed and enables flap positioning. 16Gx12 mm Cannula with two parallel extension angled tips set 1 mm, apart Each extension has four




LASIK INTERFACE CANNULA

[R] [S]




JD375-21 21G
JD375-23 23G
JD375-25 25G

Designed to wash the interface of the flap. Continuous irrigation through two side & one front holes.



LASIK IRR. CANNULA (BURRATO)

[R] [S]



JD379-23 23G
JD379-25 25G



45° Angled Triport

LASIK IRR. CANNULA (QUAD PORTED)

[R] [S]



JD381-23 23G

The Specially designed irrigation Cannula features four ports to irrigate the surgical interface and flush any debris collected during surgery. 23 GA Multi-port Irrigation tip.



LASIK FOUNTAIN CANNULA

[R] [S]



JD383-25

Both fountain cannulae are 25G(A) with a single port and a 10mm curved angel to conform to the corneal shape. The tip is a smooth closed end design slightly tapered to allow insertion under the flap. 25G(A)x1", Right 25G(A)x 1", Left


Glaucoma

LASIK CANNULA (DAYA)

[R] [S]



JD387-23 23G x 20mm



45° Angled 10 mm from end.
Flattened tip 0.30 mm

LASIK FLAP SPATULA/IRRIGATOR (HERSH)

[R] [S]



JD388-23 23G
JD388-25 25G
JD388-26 26G
JD388-27 27G

45° Angled 4 mm / 9 mm from Spatulated tip




GLAUCOMA REVISION PICK

[R]



JD391-27 27G x 16mm



60° Angle Length 7 m.m.
Suitable for Superior approach.

GLAUCOMA VALVE INTRODUCER NEEDLE

[S]



JD393

26 G/23 G Tapered tip for easy insertion.
Valve tube can be slide on 26G Needle for introduction.



VISCO CANALOSTOMY CANNULA

[R]



JD395-34 34G



Micro-Gauge Cannula for viscoelastic injection during Glaucoma Surgery. 5 mm Beveled Tip, slightly angled.

ANTERIOR CHAMBER CANNULA (MENDEZ)

[R]



JD399-27 27G x 10mm

Round tip to reform the anterior chamber
During Glaucoma Procedures Angled 45° 3 mm bend to tip




Reusable & Disposable Cannulae



Miscellaneous & New Arrivals

J301 ADOPTER




Luer lock 6% taper

J303 SILICONE BULB



J305 SILICONE TUBE WITH ADOPTER




200 mm Length

J311 MICRO SURGICAL HANDLE WITH SIDE FINGER TIP CONTROL




Standard male/female luer connection with side hole to control air with finger tip. Used with automated aspiration unit.

J317 CONNECTING TUBING W/ADOPTER (SIMCOE)




Connects 20/35 CC syringe to variety instruments

J319 IRR HANDLE WITH FINGER TIP CONTROL



Specially designed Irr. handle with flexible silicone tube.

J321 IRR. HANDLE WITH FINGER TIP CONTROL (PHACO TYPE)



Specially designed handle with silicone tube irrigation with finger tip control Phaco type.

J329 NEEDLE CORNEAL KALTS



J-329-6 6mm J-329-8 8mm



J-331 Needle Superior Rectus

Curved 3, 4, 5, 6

JD341 DISPOSABLE NEEDLE SIZE : 30G X 1/2"



Fine needle for corneal hydration, sterile. Box of 100 nos.

J335A DOUBLE MALE LUER ADAPTOR



3.5 c.m. suitable for connecting Phaco & I/A machine tubing to instruments.

IRRIGATION ASPRIATION SYSTEM

J307 MICRO SURGICAL HANDLE STAINLESS



Both sides having male luer connection

J309 MICRO SURGICAL HANDLE STAINLESS



J309 11 cm.
J309B 7.5 cm.

Standard Male/Female luer connection

J401 IRRIGATING CHOPPER



BLUNT TIP



SHARP TIP


JD205 COBRA CANNULA WITH IRRIGATING HANDLE



Bottom View




ASPIRATING CANNULA WITH POLISHER



JD357
JD357-21G (.40mm Port), JD357-23G (.30mm Port)


IRRIGATING CANNULA



JD359
JD359-21G (.50mm Port), JD359-23G (.40mm Port)

IRRIGATION ASPRIATION SYSTEM

BIMANUAL SET



TITANIUM

STEEL

DISPOSABLE

Description	Titanium	Steel	Disposable	IRRIGATION	ASPIRATION
BIMANUAL I&A SET 21G, WITHOUT POLISHER	JT451-21	J451-21	JD451-21		
BIMANUAL I&A SET 22G, WITHOUT POLISHER	JT451-22	J451-22	JD451-22		
BIMANUAL I&A SET 23G, WITHOUT POLISHER	JT451-23	J451-23	JD451-23		
BIMANUAL I&A SET 25G, WITHOUT POLISHER	JT451-25	J451-25	JD451-25		
BIMANUAL I&A SET 21G, WITH POLISHER	JT453-21	J453-21	JD453-21		
BIMANUAL I&A SET 22G, WITH POLISHER	JT453-22	J453-22	JD453-22		
BIMANUAL I&A SET 23G, WITH POLISHER	JT453-23	J453-23	JD453-23		
BIMANUAL I&A SET 25G, WITH POLISHER	JT453-25	J453-25	JD453-25		
BIMANUAL I&A SET 21G, BULLET SHAPE WITH POLISHER	JT455-21	J455-21	JD455-21		
BIMANUAL I&A SET 22G, BULLET SHAPE WITH POLISHER	JT455-22	J455-22	JD455-22		
BIMANUAL I&A SET 23G, BULLET SHAPE WITH POLISHER	JT455-23	J455-23	JD455-23		
BIMANUAL I&A SET 25G, BULLET SHAPE WITH POLISHER	JT455-25	J455-25	JD455-25		

CO-AXIAL



TITANIUM

STEEL


DISPOSABLE

Description	Titanium	Steel	Disposable
CO-AXIAL I&A HANDPIECE 19G, STRAIGHT	JT401-19	J401-19	JD401-19
CO-AXIAL I&A HANDPIECE 21G, STRAIGHT	JT401-21	J401-21	JD401-21
CO-AXIAL I&A HANDPIECE 23G, STRAIGHT	JT401-23	J401-23	JD401-23
CO-AXIAL I&A HANDPIECE 19G, ANGLED 30°	JT403-19	J403-19	JD403-19
CO-AXIAL I&A HANDPIECE 21G, ANGLED 30°	JT403-21	J403-21	JD403-21
CO-AXIAL I&A HANDPIECE 23G, ANGLED 30°	JT403-23	J403-23	JD403-23




SPECULUMS

10-100 ALFONSO SPECULUM



Infant lid speculum
with finger grips &
interlacing blades
to facilitate insertion.

10-101 BARRAQUER EYE SPECULUM



10-101 Child size
8 mm blades.
10-102 Adult size
14 mm blades
10-102A Adult size
14 mm blades
10-106 Infant size
4 mm blades

10-105 BARRAQUER SPECULUM SOLID BLADES



10-104 Child size
8 mm blades.
10-105 Adult size Temporal
14 mm blades
10-105-A Adult size, Nasal
14 mm blades

- Blade spread 18 mm
- Dull Finish
- Overall length 40 mm

10-107 KRATZ-BARRAQUER EYE SPECULUM OPEN BLADES



10-107 Adult size, Temporal
10-107 A Adult size, Nasal
10-108 Child size

- Blade spread 18 mm
- Dull Finish
- Overall length 40 mm

10-109 SIMCOE-BARRAQUER WIRE SPECULUM



- 15mm blades
- For retracting the surgical drape

SPECULUMS

10-111 MCINTYRE WIRE SPECULUM V SHAPED BLADES



- 10-111 Adult size
- Gently Curved Springs
- 'V' Shaped Blades
- Overall length : 40 mm
- 14 mm Blade

10-112 WILLIAMS SPECULUM CHILD



- 10-112 Child 11 MM BLADES, Eye speculum with locking screw
- 10-113 Adult size 15 mm Blades Eye speculum with locking screw

10-114 LESTER-BURCH SPECULUM



- 16mm Solid Blades
- Adjustable mechanism with locking nut

10-115 LANCASTER SPECULUM



- Adjustable Mechanism With Locking Screw
- 16mm Solid Blades With Flanges For Retaining Eye Lashes

10-116 SAUR EYE SPECULUM FOR INFANTS, 11MM

10-116.1 SAUR EYE SPECULUM FOR PREMATURE INFANTS, 2MM BLADES



- 10-116 11 MM BLADES
- 10-116.1 2 MM BLADES

SPECULUMS

10-117 CASTROVEIJO SPECULUM SMALL

10-118 CASTRVOEIJO SPECULUM MEDIUM

10-119 CASTRVOEIJO SPECULUM LARGE




10-117 Adjustable mechanism, 14mm blade

10-118 Adjustable mechanism, 16mm blade

10-119 Adjustable mechanism, 16 mm

10-120 GUYTON-PARK FENESTRATED SPECULUM



- Eye Speculum with Vertical suture posts and
- Canthus hook Fenestrated blades

10-125 PHACO SPECULUM

10-126 PHACO SPECULUM INFANT



10-125 Temporal

Adult size Flat branches

14mm Blade

10-128 Nasal

Adult size Flat branches

14 mm Blades

10-127 LIEBERMAN ASPIRATION SPECULUM



Six ports on each blade allow continuous aspiration for a debris - free stroma resulting in better postoperative vision
14 mm Blades

10-610 PHACO SPECULUM SOLID BLADE REVERSIBLE



Solid Blade

Adult size Flat branches

Reversible for

Temporal & Nasal Fit

14 mm Blades


SPECULUMS/RETRACTOR

10-611 CLARK EYE SPECULUM ADULT SIZE



Eye Speculum open blades
Small size

10-612 WEISS EYE SPECULUM SOLID BLADE ADULT



- Eye Speculum open blades
- Small size

10-613 NASAL SPECULUM



10-613 Adult 24 mm
10-614 Child 22 mm
10-615 Infant 20 mm

10-616 MURDOCK EYE SPECULUM



10-616
Fenestrated Blade
without locking screw
14 mm Blades

15-130 JAFFE WIRE LID RETRACTOR



15-130 Small pair
8 mm Blades
15-131 Large pair
14 mm Blades

RETRACTOR

15-135 DESMARRES LID RETRACTOR



- 15-135** Infant 7.5 mm
- 15-136** size 1 9.0 mm
- 15-137** size 2 11.0 mm
- 15-138** size 3 13.0 mm
- 15-139** size 4 15.0 mm
- Flat handle
- Polished finish
- Overall length 130 mm

15-140 AGRICOLA LACRIMAL SAC RETRACTOR, 3X3



- 15-140 Adult
- 15-140.1 Small
- Lacrimal sac retractor with 3x3 sharp prongs

15-141 STEVENSON LACRIMAL SAC RETRACTOR



- Curved 2 mm x 3 mm teeth
- Blades separated 25 mm
- Screw opening
- Overall length 90 mm

15-143 LID PLATE



- Jaeger
- Lid Plate

15-144 KNAPP LACRIMAL SAC RETRACTOR



- Blunt 4 mm x 8 mm four prong tip
- Flat serrated handle
- Polished Finish
- Overall length 130 mm

RETRACTOR


15-146 HELVESTON TISSUE RETRACTOR



- 15-146** 'Barbie' size 7 mm
- 15-146.1** 'Big Barbie' size 9 mm
- 15-146.2** 'Great Big Barbie' size 11 mm

15-148 GRAETHER COLLAR BUTTON

15-149 GRAETHER COLLAR BUTTON, ANGLED



- Round Knurled Handle
- Polished Finish
- Overall Length 130 mm

15-150 SCEPHENS ORBITAL RETRACTOR



- Flat Handle
- Wide Notch in blade
- 4.5 mm wide blade

15-151 SCHOCKET SCLERAL DEPRESSOR



- Double ended
- With pocket clip
- Overall length 130 mm

LACRIMAL & DCR INSTRUMENTS

15-152 SCHEPEN SCLERAL DEPRESSOR



- 15-152** Thimble Small
- 15-153** Thimble Medium
- 15-154** Thimble Large
- 5/7/9 mm wide tip
- Polished depressor
dull finished thimble
- Overall length 43 mm

20-160 KERRISON BONE ROUNGER FORCEPS



- 20-160** Size-1mm
- 20-161** Size-2mm
- 20-162** Size-3mm
- 20-162.1** Size-4mm
- Dull finish
- 9 mm opening

20-163 CITELLI BONE ROUNGER FORCEPS



- Bone Nibbling Ronguer
- Small

20-164 BEYER RONGEUR FORCEPS



- Single action

20-166 GASS CORNEOSCLERAL PUNCH



- 1.0 mm Dia
- Rotable Tip

LACRIMAL & DCR INSTRUMENTS

20-167 KELLY DESCOMET'S MEMBRANE PUNCH



- Produce clean 0.75 mm round hole without tissue tags
- Serrated squeeze action handle

20-168 NASAL PACKING FORCEPS



20-169 LACRIMAL SAC DOUBLE ENDED DISSECTOR



- Lacrimal Sac Dissector Sharp

20-170 WEST BONE CHISEL



- Lacrimal Sac Dissector Blunt

20-171 ARRUGA BONE TREPHINES

20-171.1 ARRUGA BONE TREPHINES W/O PIN



- Trephine't' Handle With Sharp Cutting Edge
- Cutting The Bone With Pin For Primary
- Without Pin For Secondary

LACRIMAL & DCR INSTRUMENTS

20-172 WEST BONE GAUGE



- Polished finish

20-173 FREER DOUBLE PERIOSTEAL ELEVATOR SHARP/BLUNT



- Polished finish

20-173.1 TRAQUAIR PERIOSTEAL ELEVATOR




- Polished finish

20-174 DCR MILLET



- Polished finish


20-175 BOWMAN LACRIMAL PROBE



- 20-175 Size 0000 and 000
- 20-176 size 00 and 0
- 20-177 size 1 and 2
- 20-178 size 3 and 4

LACRIMAL & DCR INSTRUMENTS


20-179 LACRIMAL PIGTAIL PROBE



Pigtail curved probes with eyes

- Round serrated handle
- Polished finish
- Overall length 140 mm



20-180 WILDER LACRIMAL DILATOR



20-180 Fine Taper
20-181 Medium Taper
20-182 Short Taper


- Blunt tip
- Round knurled handle
- Polished finish
- Overall length 100 mm

20-183 CASTROVIEJO DOUBLE ENDED LACRIMAL DILATOR SIZE 1 & 2





- Double Ended
- Fine and Medium Taper
- Overall length 130 mm

20-184 MULDOON LACRIMAL DILATOR



- Guarded Blunt Tip
- Overall length 100 mm


20-186 HARMS TRABECULOTOMY PROBE



20-186 Left
20-187 Right

LACRIMAL & DCR INSTRUMENTS

20-189 ROLLET RUGINE FOR LACRIMAL SAC




- For Lacrimal Sac

20-190 TOOKS CORNEAL KNIFE




- Corneal Knife
- Straight Blade
- Curved Cutting Edge

20-192 FRANCIS SPUD




- Foreign Body Spud Angled
- Flag Shaped Tip

20-193 "GOLF CLUB" FOREIGN BODY SPUD



- Flat Handle

25-195 GRAEFE IRIS HOOK SHARP



- Iris Hook Small


25-196 GRAEFE IRIS HOOK BLUNT



- Iris Hook Small


HOOKS

25-202 GUTHERI DOUBLE FIXATION HOOK



- Flat Handle
- Sharp Small

25-204 BURCH FIXATION PICK




- Sharp
- Double Hooks

25-207 CHAVASSE MUSCLE HOOK




- Flat handle
- Dull finish
- Overall length 140 mm

25-208 JAMESON MUSCLE HOOK



- 2 mm bulbous tip
- 6 mm flat hook
- Flat handle
- Dull finish
- Overall length 140 mm

25-210 GRAEFE MUSCLE HOOK




- 25-210 Small
- 25-211 Medium
- 25-212 Large
- Flat handle
- Dull finish
- Overall length 140 mm

25-214 GREEN STRABISMUS HOOK




- Tenotomy Hook Blunt
- Small

HOOKS/SPATULAS

25-215 STEVENS CURVED TENOTOMY HOOK



- Flat handle
- Dull finish
- Overall length 140 mm

25-217 GASS RETINAL DETACHMENT HOOK



- Flat handle dull finish
- Overall length 140 mm
- 12 mm flat hook
- 1.5 mm hole

25-218 WILLS COPPER CAUTERY (Ball Shape)
25-218.1 WILLS COPPER CAUTERY (Cone Shape)



- 12 mm copper ball
- Flat handle dull finish
- Overall length 140 mm

30-220 ELSCHNING CYCLODIALYSIS SPATULA


- Flat 1mm Wide

30-221 CASTROVIEJO DOUBLE ENDED CYCLODIALYSIS SPATULA



- 0.5 mm by 10 mm and
- 0.5 mm by 15 mm blade
- Flat handle
- Overall length 130 mm

30-222 CASTROVIEJO DOUBLE ENDED CYCLODIALYSIS SPATULA


- 0.25 mm
- Flat handle
- Overall length 130 mm


SPATULAS

30-223 DASTOOR PUPIL & IRIS REPOSITOR




- Flat handle
- Overall length 145 mm

30-225 KNAPP SPATULA & HOOK




- Double-ended Spatula

30-226 MCPHERSON IRIS SPATULA



- 1mm wide
- 2mm wide
- Overall length 120 mm

30-227 WHEELER DOUBLE-ENDED SPATULA LIGHTY




- Double-ended Spatula Lightly Curved

30-228 GIRALD SYNECHIAE SPATULA SEMISHARP TIP



- 0.5mm Wide With Semi-sharp Tip

30-229 OBSTBAUM LENS SPATULA ANGLED BLADE



- Lens Spatula Angled Blade
- 0.25mm Wide

30-230 JAFFE LENS SPATULA



- Width 0.50 mm, length 12 mm
- Flat handle
- Overall length 125 mm

LOOPS

- 30-231 KOCH NUCLEUS SPATULA**
30-232 KOCH NUCLEUS SPATULA WITH NOTCHES



- 30 mm long narrow spatula blade
- Round handle
- Overall length 120 mm

- 30-233 PATON DOUBLE ENDED SPATULA**




- Paton
Spatula & Spoon
For Corneal Button
Transfer
- Overall length 130 mm

- 30-234 GILL-WELSH LENS LOOP & OLIVE TIP PROBE**




Lens Loop And Probe
With Serrations

- 35-239 NEW ORLEANS EYE & EAR LENSE LOOP**




- Lens Loop And Probe
With Serrations

- 35-241 WEBER LENS LOOP**




- 35-241 Small**
35-242 Large
- Rectangular handle
 - Overall length 130 mm

- 35-243 GILLS-WELSH THIN BLADE LENSE LOOP**



- Lens Loop And Probe


- 35-245 ARIT LENSE LOOP**



- Lens Loop
- 3mm Wide Serrated Strong
- Curve for Small Incision Nucleus
Fragment Removal


LOOPS/ SPOONS & CURETTES

35-247 LENS LOOP WITH SERRATION



- Rectangular handle
- Overall length 130 mm

40-250 BUNGE EVISCERATION SPOON



40-250 size 1 Small
40-251 size 2 Large

- Flat handle
- Overall length 130 mm

40-252 WELLS ENUCLEATION SPOON




- Flat handle
- Overall length 145 mm

40-254 BARRAQUER LENTICULAR SPOON



- Flat handle
- Overall length 145 mm

40-255 MEYERHOEFFER CHALAZION CURETTE



- 40-255 size 0 1.5 mm**
- 40-256 size 1 1.75 mm**
- 40-257 size 2 2.0 mm**
- 40-258 size 3 2.50 mm**
- 40-259 size 4 3.50 mm**
- Flat handle
- Overall length 130 mm

40-261 NIGHTINGALE CAPSULE POLISHER




- Capsule Polisher
- 1.75mm dia. Open Ring
- Semi-Sharp Edge


NUCLEUS ROTATORS/MANIPULATORS

45-270 BECHERT NUCLEUS ROTATOR



- "Y" Shape tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

45-271 JAFFE-BECHERT NUCLEUS ROTATOR



- Angled Lens Pusher With Vertical Y-shaped Tip

45-273 GIRALD SYNECHIAE SPATULA


- Synechiae Spatula
Semi-sharp Tip 0.5mm
Wide, Angled

45-274 LIPPMAN SIDEPORT NUCLEUS SPATULA


- 0.5mm Diameter Round Shaft
With Blunt Forked Tip


45-275 HIRSCHMAN NOTCHED LENS SPATULA


- Lens Spatula Flat
• 0.5mm Wide With Small
Notches For Manipulating IOL

45-278 LESTER LENS MANIPULATOR

**45-277 Straight
45-278 Angled**

- Round Straight shaft
with 0.2 mm bulbous tip
- Round serrated handle
- Dull finish
- Overall length 130 mm


45-280 AKER LENS PUSHER STRAIGHT


- Lens Pusher With Smooth
Fish Mouth Shaped Tip

NUCLEUS ROTATORS/MANIPULATORS




45-282 BURRATO HOOK




- Bulbous tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

45-283 BURRAQUER HOOK FOR ICL 45-283.1 ICL MANIPULATOR




- Flat tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

45-287 SINSKEY HOOK




- 10 mm from tip to bend
- Round serrated handle
- Dull finish
- Overall length 130 mm

45-295 KUGLEN HOOK



- 45-294 Straight
- 45-295 Angled
- "H" shape tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

45-296 KUGLEN-DE-LA-VEGA LOL HOOK IOL MANIPULATOR




- Clover Leaf Model Iris Hook
Lens Manipulator

NUCLEUS ROTATORS/MANIPULATORS




45-297 HIRSCHMAN IRIS/LENS HOOK ANGLED




- Lens Spatula
- Flat, 0.5mm Wide
- Small notches for Manipulating IOLs

45-299 FENZL HOOK




- 45-298 Straight
- 45-299 Angled
- Overall length 130 mm

45-300 MALTZMAN-FENZL LENS MANIPULATOR HOOK



- Modified No Hole Lens Manipulator
- 0.15mm diameter tip

45-303 MALONEY LENS MANIPULATOR STRAIGHT




- IOL Manipulator Vaulted Shank 0.25mm Diameter Tip

45-304 LITTLE LENS MANIPULATOR



- No Hole Lens Manipulator tip
- Width 1.20mm STRAIGHT

45-305 NUCLEUS SUSTAINER




- Ball tip
- Round serrated handle
- Dull finish
- Overall length 130 mm

NUCLEUS ROTATORS/MANIPULATORS




45-306 LEWICKY HOOK



- 0.2 mm hook on curved 10 mm shaft
- Round serrated handle
- Dull finish
- Overall length 130 mm

45-307 PISACANO NUCLEUS ROTATOR



- Nucleus Rotator With Flat Disc Shaped Tip

45-308 TENNANT NUCLEUS BALL TIP MANIPULATOR




- Button shape Tip
- Dull finish
- Round knurled handle
- Overall length 130 mm

45-311 KANSAS NUCLEUS TRISECTOR




- Round knurlend handle
- Dull finish
- Overall length 130 mm

45-312 KANSAS NUCLEUS BISECTOR



- Bottom edged shaft
- Round knurled handle
- Dull finish
- Overall length 130 mm

45-313 KANSAS NUCLEUS VECTIS SOLID




- Round serrated handle
- Dull finish
- Overall length 115 mm

CHOPPER


45-314 PHACO CHOPPER



45-314 Sharp



45-314.1 Blunt



45-314.2 Nagahara




45-314.3 Agarwal




45-314.4 Verghese



45-314.8 Mohan Rajan



45-314.9 Arul




45-314.5 Olson




SPATULA/FIXATION RINGS

45-323 FUKASAKU LASIK SPATULA



- Lasik Spatula for Lifting & Dissecting the Corner of Flap
- Fixed round handle
- Overall length 135 mm

45-324 GUIMARAES LASIK SPATULA



- Malleable
- Fixed round handle
- Overall length 135 mm

45-325 GLOBE FIXATION RING FOR PRK/LASIK PROCEDURE



- Fixed handle
- Inside Dia. 12 mm
- Overall length 135 mm

45-326 GLOBE FIXATION RING




- open universal
- 12 mm central opening
- 3/4 open ring
- Fixed round handle
- Overall length 135 mm

45-327 OSHER-NEUMANN RADIAL MARKER



- Available with 4-12 blades
- Fixed round handle
- Overall length 135 mm

45-328 HOFFER OPTIC ZONE MARKER CROSS HARIS



- Diametre 2.5 mm to 10 mm
- Fixed round handle
- Overall length 135 mm

FIXATION RINGS

45-329 THORNTON OPTIC ZONE MARKER WITH CENTRE



- Size 3mm - 9 mm
- Round Handle

45-330 GIMBEL "C" STYLE FIXATION RING




45-330 16 mm diameter ring with teeth

45-331 13 mm diameter ring with teeth

- 3/4 open ring with swivel
- Round Handle
- Overall length 135 mm

45-682 FLIERINGA FIXATION RING



- Scleral Fixation Ring
- Set of 8 Sizes (15-22mm)

50-335 FLEXIBLE STAINLESS STEEL RULER



- Stainless Steel

CALIPERS, GAUGES

50-336 CASTROVIEJO MARKING CALIPER, STRAIGHT

50-337 CASTROVIEJO MARKING CALIPER, CURVED



- Measures from 0 to 20 mm
- Scale reads on both sides
- Polished finish
- Polished finish
- Overall length 87 mm

CALIPERS, GAUGES


50-338 MENDEZ DEGREE GAUGE



- For accurate axial alignment of corneal incision to correct astigmatism 12mm I.D.

50-339 FIXED CALLIPER

50-339.1 FIXED CALIPER 3.0-3.5MM




- Size 3.5 mm to 4 mm
- Marks the distance from the limbus to sclera

50-806 YEOH TORIC GAUGE II 12MM I.D./14MM O.D.




- Degree Gauge For Accurate Axial Alignment Of Corneal Incision To Correct Astigmatism 12mm I.d

50-808 GREATHER HALF CIRCLE AXIS MARKER




FORCEPS

55-354 BISHOP-HARMON DRESSING FORCEPS, SERRATED




- Criss Cross Serrated tips
- Flat handle
- Dull finish
- Overall length 90 mm

55-355 BONN FORCEPS, SHORT HANDLE




- 1x2 teeth, 0.12 mm
- Flat short handle
- Dull finish
- Overall length 75 mm

55-356 BONN FORCEPS, 1X2 12MM TEETH, TYING





- Very delicate
1x2 teeth
0.12mm

55-357 BONN FORCEPS LONG HANDLE





- Flat long Serrated handle
- 1x2 teeth, 0.12 mm
- Dull finish
- Overall length 110 mm

55-358 HESS IRIS FORCEPS, SWISS MODEL VERY DEL.



- Iris Forceps Swiss Model
1 X 2 Teeth


55-359 HESS IRIS FORCEPS, STANDARD MODEL VERY



- Iris Forceps Standard Model
1 X 2 Teeth


FORCEPS

55-362 WILLS HOSPITAL UTILITY FORCEPS, SERRATED



- Utility Forceps
- Criss Cross Serrations

55-363 BONACCOLTO UTILITY FORCEPS, SERRATED




- Flat Serrated Handle
- Long Serrations

55-364 GILL IRIS FORCEPS




- Curved Shafts
- Criss Cross Serrations
- Flat handle

55-365 GILL-BOTVIN IRIS FORCEPS




- Iris Forceps Curved
- Criss-cross Serrated Tips

55-370 ELSCHING FIXATION FORCEPS 5MM JAW



- Suturing Forceps
- 1 X 2 Teeth
- Tying Platform
- 0.12mm

55-371 GREEN FIXATION FORCEPS 5MM JAW



- Fixation Forceps
- Multiple Teeth
- 5mm Wide Jaws

FORCEPS

55-340 SERREFINE



55-340 Straight small size

55-341 Curved small size

55-342 HALSTED HEMOSTATIC MOSQUITO FORCEPS




55-342

- Straight serrated jaws
- Ring handle
- Dull finish
- Overall length 120 mm

55-343

- Curved serrated jaws
- Ring handle
- Dull finish
- Overall length 120 mm

55-344 HARTMAN HEMOSTATIC MOSQUITO FORCEPS



55-344 Straight

- Straight serrated jaws
- Ring handle
- Dull finish
- Overall length 75 mm

55-345 Curved

- Curved serrated jaws
- Ring handle
- Dull finish
- Overall length 75 mm

55-346 TOWEL CLAMP - CROSS ACTION



- Flat handle
- Polished finish
- Overall length 85 mm


55-347 BACKHAUS TOWEL CLAMP



- Towel Clamp
- 3.5" (90mm) Long
- 5.25" (135mm) Long


FORCEPS

55-348 STRAIGHT DRESSING FORCEPS WITH DELICATE SERRATIONS




- Straight 6 mm serrated tips
- Dull finish
- Flat serrated handle
- Overall length 100 mm

55-349 CURVED DRESSING FORCEPS WITH SERRATIONS




- Curved serrated tips
- Flat serrated handle
- Dull finish
- Overall length 100 mm

55-350 IRIS FORCEPS WITH TEETH




- Straight shaft
- 1x2 teeth
- Flat serrated handle
- Dull finish
- Overall length 100 mm

55-351 IRIS FORCEPS WITH TEETH




- Curved shaft
- 1x2 teeth
- Dull finish
- Overall length 90 mm

55-352 SWISS MODEL IRIS FORCEPS, CURVED VERY



- Model Iris Forceps Curved
Very Delicate 1 X 2 Teeth


55-353 BISHOP-HARMON FORCEPS



- 3 mm tying platform
- 1x2 teeth 0.3 mm
- Flat handle
- Dull finish
- Overall length 90 mm


FORCEPS

55-372 GREEN FIXATION FORCEPS 10MM JAW




- Fixation Forceps
- Multiple Teeth
- 10mm Wide Jaws

55-373 FIXATION FORCEPS




- 1x2 teeth 0.5 mm
- Straight shafts
- Serrated handle
- Dull finish
- Overall length 115 mm

55-374 CASTROVIEJO FIXATION FORCEPS, 2X3 TEETH




- Fixation Forceps Multiple Teeth
- 5mm Wide Jaws

55-375 LESTER FIXATION FORCEPS




- 1x2 Teeth 0.5 mm
- Straight shafts
- Serrated Short handle
- Dull finish
- Overall length 90 mm

55-376 LESTER FIXATION FORCEPS, 2X3 TEETH 9CM



- Straight, 2x3 Teeth 9 Cm


55-377 BORES U-SHAPED FIXATION FORCEPS



- Fixation Forceps 3mm Spread Each 0.12mm
- 1 X 2 Teeth Straight

FORCEPS


55-378 BORES U-SHAPED FIXATION FORCEPS, ANGLED



- Fixation Forceps 3mm Spread Each 0.12mm
- 1 X 2 Teeth Curved

55-379 ARROWSMITH FIXATION FORCEPS, STRAIGHT


55-380 ARROWSMITH FIXATION FORCEPS, ANGLED



- Corneal Fixation Forceps With Lock "u" Shaped
- 13mm Spread, Each 0.12mm
- 1 X 2 Teeth Straight With Catch

55-381 KREMER FIXATION FORCEPS

55-382 KREMER ANGLED FIXATION FORCEPS



- Tips separated by 13 mm
- 0.12 mm, 1x2 teeth
- Flat handle
- Dull finish
- Overall length 105 mm


55-383 DASTOOR SUPERIOR RECTUS FORCEPS



- 1x2 Teeth
- Serrated handle
- Dull finish, Overall length 105 mm


FORCEPS

55-385 JEWELLER FORCEPS




- Straight shafts
- Pointed tips
- Tying Platform
- Smooth handle
- Dull finish
- Overall length 110 mm

55-386 JEWELLER FORCEPS




- Straight shafts
- Pointed tips
- Tying Platform
- Smooth handle
- Dull finish
- Overall length 110 mm

55-387 JEWELLER FORCEPS




- Straight shafts
- Pointed tips
- Tying Platform
- Smooth handle
- Dull finish
- Overall length 110 mm

55-388 JEWELLER FORCEPS




- Straight shafts
- Pointed tips
- Tying Platform
- Smooth handle
- Dull finish
- Overall length 110 mm

55-389 JEWELLER FORCEPS



- Straight shafts
- Pointed tips
- Tying Platform
- Smooth handle
- Dull finish
- Overall length 110 mm


55-392 JAMESON MUSCLE FORCEPS - RIGHT



- 12 mm angled shafts
- Scaled Jaws
- Teeth 1 mm on one jaw
opp. holes on the other jaw
- Serrated handle
- Dull finish
- Overall length 110 mm

FORCEPS

55-393 JAMESON MUSCLE FORCEPS - LEFT



- 12 mm angled shafts
- Scaled Jaws
- Teeth 1 mm on one jaw opp. holes on the other jaw
- Serrated handle
- Dull finish
- Overall length 110 mm

55-394 BERKE PTOSIS FORCEPS



55-394 20 mm
55-395 27 mm

- With side lock
- Longitudinally grooved jaws


55-396 SNELLON ENTROPIUM FORCEPS



55-396 Right
55-397 Left

- Solid lower plate, open upper plate
- Locking thumb screw mechanism
- Polished finish
- Overall length 90 mm

55-398 LAMBERT CHALAZION FORCEPS SMALL




- Round solid lower plate open upper plate
- 9 mm inside diameter of upper plate
- Locking thumb screw mechanism
- Dull finish
- Overall length 97 mm

55-399 LAMBERT CHALAZION FORCEPS MEDIUM



- Round solid lower plate open upper plate
- 12 mm inside diameter of upper plate
- Locking thumb screw mechanism
- Dull finish
- Overall length 97 mm

55-400 DESMARRES CHALAZION FORCEPS



55-400 Small 20mm
55-401 Medium 26mm
55-402 Large 31mm

- Round solid lower plate open upper plate
- 20/26/31 inside diameter of upper plate
- Locking thumb screw mechanism
- Dull finish
- Overall length 97 mm


FORCEPS

55-404 ERHARDT LID FORCEPS



- Lid Forceps

55-405 HUNT CHALAZION FORCEPS




- Chalazion Forceps Round
- 12mm Diameter

55-406 WEIS CHALAZION FORCEPS




- Grooved Upper Jaw

55-408 BARRAQUER CILIA FORCEPS



- Curved Tips
- Smooth jaws
- Flat handle
- Dull finish
- Overall length 100 mm


55-409 BEER CILIA FORCEPS



- Smooth jaws
- Serrated handle
- Dull finish
- Overall length 80 mm


FORCEPS

55-410 ZIEGLER FORCEPS CILIA




- Smooth jaws
- Serrated handle
- Dull finish
- Overall length 94 mm

55-411 HUNT ARRUGA CAPSULE FORCEPS




- Capsule Forceps

55-412 HUNT ARRUGA CAPSULE FORCEPS NEW CURVE




- Capsule Forceps New Curved

55-420 JAFFE TYING FORCEPS STRAIGHT




- For 8-0 to 11-0 sutures
- 6 mm tying platform
- Flat handle
- Overall length 100 mm

55-421 JAFFE TYING FORCEPS CURVED



- For 8-0 to 11-0 sutures
- 6 mm tying platform
- Flat handle
- Overall length 100 mm


55-422 SINSKY TYPING FORCEPS STRAIGHT



- Tying Forceps Very Delicate Straight
- Smooth Jaws Ideal For 8-0 To 11-0 Sutures


FORCEPS

55-423 SINISKY TYING FORCEPS CURVED




- Tying Forceps Very Delicate Curved
- Smooth Jaws Ideal For 8-0 To 11-0 Sutures

55-424 TENNANT TYING FORCEPS STRAIGHT




- Extra delicate smooth jaws
- 6 mm tying platform
- 6 mm For 9-0 to 11-0 sutures
- Round serrated handle
- Overall length 100 mm

55-425 TENNANT CURVED TYING FORCEPS




- Extra delicate smooth jaws
- 6 mm tying platform
- 6 mm For 9-0 to 11-0 sutures
- Round serrated handle
- Overall length 100 mm

55-429 MCPHERSON STRAIGHT FORCEPS




- Smooth jaws
- 5 mm tying platform
- For 8-0 to 11-0 sutures
- Flat handle
- Overall length 85 mm

55-430 MCPHERSON STRAIGHT FORCEPS LONG HANDLE



- Straight shafts
- Pointed tips
- Tying platform
- Polished finish
- Overall length 110 mm


55-431 MCPHERSON TYING FORCEPS ANGLED



- Angled shafts
- Pointed tips
- Tying platform 5 mm
- Serrated handle
- Dull finish
- Overall length 85 mm


FORCEPS

55-432 MCPHERSON TYING FORCEPS LONG HANDLE ANGLED




- Straight shafts
- Pointed tips
- Tying platform 5mm
- Flat handle
- Polished finish
- Overall length 110 mm

55-433 HARMS TYING FORCEPS STRAIGHT




- For 7-0 to 8-0 sutures
- 5 mm tying platform
- Straight Shafts
- Pointed tips
- Serrated handle
- Overall length 100 mm

55-434 HARMS TYING FORCEPS CURVED




- For 7-0 to 8-0 sutures
- 5 mm tying platform
- Curved Shafts
- Pointed tips
- Serrated handle
- Overall length 110 mm

55-435 CASTROVEIJO SUTURE FORCEPS




- 55-435 1x2 Teeth 0.12 mm
- 55-436 1x2 Teeth 0.20 mm
- 55-437 1x2 Teeth 0.30 mm
- 55-438 1x2 Teeth 0.50 mm
- 55-439 1x2 Teeth 0.90 mm
- 55-440 1x2 Teeth 1.50 mm
- With tying platform
- Overall length 110 mm

55-441 CASTROVEIJO ANGLED SUTURING FORCEPS



- 1x2 teeth 0.12 mm
- Angled Shaft
- 6 mm tying platform
- Dull Finish
- Flat handle
- Overall length 110 mm


55-442 BARRAQUER COLIBRI SUTURE FORCEPS W/O TYING PLATFORM



- 1x2 teeth 0.12 mm
- Flat serrated handle
- Without Tying platform
- Dull Finish
- Overall length 90 mm


FORCEPS

55-443 BARRAQUER-COLIBRI SUTURE FORCEPS, 1X2 TEETH, WITH TYPING PLATFORM




- Suture Forceps
- 1x2 Teeth With Typing Platform

55-444 HARMS-COLIBRI SUTURE FORCEPS




- 1x2 teeth 0,12 mm
- Flat handle
- Tying platform
- Dull Finish
- Overall length 90 mm

55-445 GILLS COLIBRI CAPSULE FORCEPS SMOOTH




- Capsule Forceps
- Very Fine
- Pointed Tips

55-446 PIERSE COLIBRI FIXATION FORCEPS




- Cup Shaped tips
- Tying platforms
- Dull Finish
- Flat serrated handle
- Overall length 90 mm

55-448 CASTROVEIJO-COLIBRI FORCEPS



- 1x2 Teeth, 0.12 mm
- Fine Delicate
- Flat handle


55-449 TENNANT-COLIBRI CROSS-ACTION FORCES



- Colibri Forceps With Tying Platform
- 1 X 2 Teeth, 0.12mm


FORCEPS

55-452 NUGENT UTILITY FORCEPS




- Smooth Jaws

55-453 PLAUFIQUE SUTURING FORCEPS




- 6 mm tying platforms
- Flat serrated handle
- 0.12 mm, 1x2 Teeth
- Dull Finish
- Overall length 87 mm

55-454 ALABAMA FORCEPS




- Tying Forceps Very Delicate
- Smooth Jaws Ideal For 8-0 To 11-0 Sutures Curved

55-455 MOORFIELD SUTURING FORCEPS




- Conjunctival suture forceps
- Horizontally serrated jaws
- With groove on tip
- Flat serrated handle

55-456 LIM'S SUTURING FORCEPS



- 1x2 extra fine teeth
- Tying platforms
- Dull Finish
- Flat serrated handle
- Overall length 87 mm


55-457 MCPHERSON SUTURING FORCEPS



- 1x2 extra fine teeth
- Tying platforms
- Dull Finish
- Flat serrated handle
- Overall length 90 mm


FORCEPS

55-458 MCPHERSON SUTURING CORNEAL FORCEPS 1X2 FINE TEETH ANGLED




- Corneal Forceps 7.5mm Long Angled Jaws
- 1 X 2 Teeth, 0.3mm

55-459 HOFFER-MCPHERSON FORCEPS 11MM ANGLED WITH 0.12MM 1X2 TEETH




- Corneal Forceps Very Delicate 11mm Long Angled Jaws
- 1 X 2 Teeth, 0.12mm

55-460 THORPE CORNEAL FORCEPS, 2X3 0.12MM




- 2x3 Teeth, 0.12mm With
- Tying Platform

55-461 GIRARD CORNEAL FORCEPS 1X2 TEETH 0.12




- 1x2 Teeth, 0.12mm With
- Tying Platform

55-462 KELMAN-MCPHERSON FORCEPS



- Angled shaft
- 8mm Tying platforms
- Dull Finish
- Flat handle
- Overall length 100 mm


55-463 BECHERT MCPHERSON FORCEPS



- Angled shaft
- 10mm Tying platforms
- Dull Finish
- Flat handle
- Overall length 100 mm


FORCEPS

55-464 GASKIN FRAGMENT FORCEPS, 11MM SMOOTH




- 11mm long smooth jaws
- Angled,

55-468 BECHERT IOL IMPLANTATION FORCEPS




- Curved shaft
- Flat handle
- Dull Finish
- Overall length 120 mm

55-469 SHEPHARD IOL FORCEPS




- Lens Holding forceps
- 3 mm lower & 1.5 mm Upper jaws
- Curved shaft
- Flat handle
- Dull Finish
- Overall length 120 mm

55-470 CLAYMAN LENS HOLDING FORCEPS




- Curved shaft
- Angled 3 mm tips
- Flat handle
- Dull Finish
- Overall length 120 mm

55-471 BLADYES IOL INSERTION FORCEPS



- Angled shaft
- Flat handle
- Dull Finish
- Overall length 120 mm


55-475 DODICK LENS HOLDING FORCEPS, ANGLED TO LEFT



- Lens Holding Forceps
- Angled Smooth Jaws
- Angled Forward For Comfortable Hand Position

FORCEPS

55-476 DODICK LENS HOLDING FORCEPS, ANGLED TO RIGHT




- Lens Holding Forceps
- Angled Smooth Jaws
- Angled Forward For Comfortable Hand Position

55-478 ANIS LENS HOLDING FORCEPS




- Angled shaft
- Round handle
- Dull Finish
- Overall length 120 mm

55-479 SINSKEY LENS HOLDING FORCEPS




- Straight shaft
- Delicate T-shaped Lower Jaw
- Flat handle
- Dull Finish
- Overall length 120 mm

55-480 INAMURA CAPSULORHEXIS FORCEPS




55-481 UTRATA CAPSULORHEXIS FORCEPS



- 11mm angled tip
- With iris stop platform
- Flat handle
- Overall length 100 mm


55-481.1 UTRATA CAPSULORHEXIS FORCEPS, STRAIGHT JAWS, SMALL FLAT HANDLE



- Straight Jaws
- Small Flat Handle


FORCEPS

55-482 UTRATA CAPSULORHEXIS FORCEPS STRAIGHT JAWS ROUND HANDLE




- 11mm angled tip
- With iris stop platform
- Round handle
- Overall length 100 mm

55-483 MASKET CAPSULORHEXIS FORCEPS CURVED JAWS FLAT HANDLE




- 11mm curved tip
- With iris stop platform
- Flat handle
- Overall length 100 mm

55-484 DODICK NUCLEUS CRACKER




- Cross action
- Round serrated handle
- Overall length 100 mm

55-485 KANSAS NUCLEUS REMOVAL FORCEPS



- Angled shaft with fine teeth
- For extracting nucleus fragments in SICS
- Flat handle

55-486 FAULKNER LENS HOLDING FORCEPS






- Smooth round jaws
- For Used to safely pick foldable IOL from case
- Flat handle

55-487 FAULKNER LENS FOLDING FORCEPS




- Lens Inserting Forceps
For Folding And Inserting
Soft IOLs Through Small Incision


FORCEPS

55-489 ERNEST NUCLEUS CRACKER

55-490 MC DONALD SOFT IOL INSERTING FORCEPS CROSS-ACTION

55-491 SPALEEK SOFT IOL INSERTION FORCEPS



- For inserting IOL through small incision
- Flat handle
- Overall length 100 mm

55-493 PIERSE TYPE MICRO FORCEPS


- Pierse Type
- Micro Forceps
- Straight / Curved
- Long Handle

55-497 FOLDING FORCEPS FOR ACRYLIC & SILICONE



- Angled shaft
- Flat handle
- Dull finish
- Overall length 100 mm

55-498 UNIVERSAL PREFOLDER FOR & SILICONE


- Angled, cross action shaft
- Round handle
- Dull finish
- Overall length 100 mm


FORCEPS

55-619 BURATTO LASIK FLAP HOLDING FORCEPS




- Serrated & abrasive blasted inner tip
- Flat handle

55-622 PIERSE HOSKIN FORCEPS




- 55-622** 0.3 mm
- 55-623** 0.8 mm
- 55-624** 1.5 mm

55-625 AKASHOSHI PHACO PRECHOPPER STRAIGHT




- Cross action
- Round handle
- Maximum opening 2.8 mm
- Overall length 100 mm

55-626 AKASHOSHI PHACO PRECHOPPER ANGLED




- Cross action
- Round handle
- Maximum opening 2.8 mm
- Overall length 100 mm

55-627 AKASHOSHI PHACO PRECHOPPER CURVED



- Cross action
- Round handle
- Maximum opening 2.8 mm
- Overall length 100 mm



55-628 AKASHOSHI PHACO PRECHOPPER COMBO




- Cross action
- Bottom cutting edge on tip
- Round handle
- Maximum opening 2.8 mm
- Overall length 100 mm

FORCEPS


55-656 SCLERAL PLUG HOLDING FORCEPS



55-657 MOODY FIXATION FORCEPS RIGHT



55-658 MOODY FIXATION FORCEPS LEFT




- 1x2 teeth, 0.5 mm with lock
- Overall length 100 mm

55-659 IRIS CLAW LENS FORCEPS WITHOUT LOCK



55-659.1
Iris Claw Lens Forceps
With Lock

55-681 NEVYAS CAPSULORHEXIS FORCEPS, EYSTOTOME SHAPED TIPS, EXTREMELY THIN



Nevvas
Capsulorhexis Forceps Extremely
Thin 11mm Long Shanks
Cystotome Shaped Tips For
Penetrating And Grasping
Anterior Capsule

55-720 MEIBOMIAN GLAND DEPRESSING FORCEPS




FORCEPS

55-723 MASTROTA MEBOMIAN GLAND FORCEPS




55-724 GROOVED FORCEPS FOR PUNCTUM PLUGS




55-725 SPONGE HOLDING FORCEPS




55-726 DMEK GRAFT CROSS ACTION FORCEPS



55-727 MEMBRANE PEELING FORCEPS STRAIGHT




55-728 REVERSE CAPSULORHEXIS FORCEPS




FORCEPS

55-729 BIPOLAR FORCEPS, REGULAR




- Straight

55-730 BIPOLAR FORCEPS, MCPHERSON ANGLED




- Angled

55-731 BIPOLAR FORCEPS, MCPHERSON STRAIGHT




- Straight


55-821 SERRATED FORCEPS CURVED SHAFT, VERTICAL ACTION



55-823 ENCLAVATION FORCEPS




80-816 GLUED IOL FORCEP




FORCEPS

55-800 SQUEEZE HANDLE FORCEP DSAEK




55-831 ICL LOADING FORCEP



VITREORETINAL FORCEPS




Straight Smooth Tip




55-494

Foreign Body




55-495

Vertical Swiss Action




55-496

Capsulorhexis




55-499

1x2 teeth, 0.12 mm




55-617

Straight Serrated




55-618

ILM




55-631

PICK



55-632

End Gripping





55-633

- Squeeze Action Handle
- 21G / 23G / 25G

SCISSORS


60-500 STRAIGHT IRIS SCISSORS



60-500 9 cm straight
60-501 9 cm curved
60-502 11.5 cm curved


- Pointed tips 28 mm
- Ring handle
- Dull finish
- Overall length 90 mm

60-503 KNAPP STRABISMUS SCISSORS



60-503 Straight
60-504 Curved

60-505 STEVENS TENOTOMY SCISSORS



60-505 Straight
60-506 Curved

- Short Model
- Blunt
- Ring handle
- Dull finish

SCISSORS

60-507 SPENCER STICH SCISSORS




- Short Model
- Blunt
- Ring handle
- Dull finish

60-508 STICH SCISSORS, RIBBON TYPE, NEEDLE POINTS SLIGHT CURVED



- Needle Points
- Light Curved
- Ring Handle

60-510 STEVENS TENOTOMY SCISSORS, RIBBON TYPE



- Tenotomy Scissors
- Blunt Tips
- Ring Handle

60-512 ENUCLEATION SCISSORS, MEDIUM CURVE



- Medium Curve
- Round blades 45 mm
- Blunt tips
- Ring handle
- Dull finish
- Overall length 148 mm


60-513 ENUCLEATION SCISSORS, STRONG CURVE



- Strong Curve
- Round blades 48 mm
- Blunt tips
- Ring handle
- Dull finish
- Overall length 148 mm


SCISSORS

60-514 BARRAQUER IRIS SCISSORS




- Blunt tips
- 7 mm blades
- Flat serrated squeeze action handle
- Dull finish
- Overall length 54 mm

60-515 VANNAS SCISSORS, STRAIGHT




- 6 mm blades
- Sharp pointed tips
- Dull finish
- Overall length 90 mm

60-516 VANNAS SCISSORS, CURVED




- 6 mm blades
- Sharp pointed tips
- Dull finish
- Overall length 90 mm

60-517 VANNAS SCISSORS, ANGLED ON SIDE




- 6 mm blades
- Sharp pointed tips
- Dull finish
- Overall length 85 mm

60-518 GILLS-VANNAS SCISSORS, ANGLED



- 6mm blade
- Sharp pointed tips
- Dull finish
- Flat handle
- Overall length 85 mm


60-519 GILLS-VANNAS SCISSORS, ANGLED



- 11mm blade
- Sharp pointed tips
- Dull finish
- Flat handle
- Overall length 85 mm


SCISSORS

60-520 MCPHERSON-VANNAS SCISSORS, CURVED




- 11 mm blades
- Flat handle
- Dull finish
- Overall length 84 mm

60-521 WESTCOTT TENOTOMY STITCH SCISSORS




- 11 mm blades
- Blunt tips
- Flat handle
- Dull finish
- Overall length 115 mm

60-522 WESTCOTT STITCH SCISSORS




- 16 mm blades
- Sharp pointed tips
- Dull finish
- Overall length 115 mm

60-523 WESTCOTT CURVED TENOTOMY SCISSORS



- 60-523 Right**
60-524 Left
- Blunt tips
 - Standard blades
 - Dull finish
 - Overall length 115 mm

60-525 MCPHERSON-WESTCOTT STITCH SCISSORS BLUNT TIPS, SMALL BLADES



- Blunt tips
- 9 mm blades
- Flat serrated handle
- Dull finish
- Overall length 105 mm

SCISSORS

60-526 JAFFE STITCH SCISSORS, MEDIUM BLADES



- Jaffe
- Stitch Scissors
- Very Sharp Pointed Tips
- Medium Blades

60-527 CASTROVEIJO CORNEAL SCISSORS, SMALL BLADES



- Curved on Flat
- Small blades 7 mm
- Flat Serrated handle
- Dull finish
- Overall length 100 mm

60-528 CASTROVEIJO CORNEAL SCISSORS, MEDIUM BLADES




- Curved on Flat
- Medium blades 9 mm
- Flat Serrated handle
- Dull finish
- Overall length 100 mm

60-529 CASTROVEIJO CORNEAL SCISSORS, LARGE BLADES



- Curved on Flat
- Large blades 11 mm
- Flat Serrated handle
- Dull finish
- Overall length 100 mm

60-530 MICRO SCISSORS STRAIGHT



- Universal miniature
- mode with blunt tips
- 11 mm blades
- Flat handle
- Dull finish
- Overall length 100 mm

60-532 CASTROVEIJO CORENEAL SCISSORS BLUNT TIP



- 9mm blades
- Flat handle
- Blunt tips
- Dull finish
- Overall length 100 mm

SCISSORS

60-533 CASTROVEIJO CORENEAL SCISSORS BLUNT TIP



60-533 11 mm blades
medium curved

60-534


- Strong curve 11 mm blades
- Flat handle
- Blunt tips
- Dull finish
- Overall length 100 mm

60-535 CASTROVEIJO CORNEAL SCISSORS, SCISSORS RIGHT




- Lower blade 0.5 mm longer than upper blades
- 9 mm blades
- Flat handle
- Dull finish
- Overall length 106mm

60-536 CASTROVEIJO CORNEAL SECTION, SCISSORS LEFT



- Lower blade 0.5 mm longer than upper blades
- 9 mm blades
- Flat handle
- Dull finish
- Overall length 106mm

60-537 CASTROVEIJO CURVED CORNEAL SCISSORS RIGHT



- Right
- 11 mm blades
- Blunt tips
- Flat handle
- Dull finish
- Overall length 106 mm


60-538 CASTROVEIJO CURVED CORNEAL SCISSORS LEFT



- Left
- 11 mm blades
- Blunt tips
- Flat handle
- Dull finish
- Overall length 106 mm

SCISSORS

60-543 TROUTMAN-KATZIN CORNEAL SCISSORS RIGHT




- Small blades
- With Stop Lock
- Overall length 100mm

60-544 TROUTMAN-KATZIN CORNEAL SCISSORS, LEFT



- Small blades
- With Stop Lock
- Overall length 100mm

60-545 KATZIN CURVED CORNEAL TRANSPLANT



- Small blades
- With Stop Lock
- Overall length 100mm

60-546 KATZIN CURVED CORNEAL TRANSPLANT



- Small blades
- With Stop Lock
- Overall length 100mm


60-547 CASTROVEIJO KERATOPLASTY SCISSORS, SIDE ANGLED



- Universal with blunt angled blades
- Lower blade 1 mm longer than upper blades
- 16 mm blades
- Flat handle
- Dull finish
- Overall length 110 mm

SCISSORS


60-548 CAPSULOTOMY SCISSORS



- Universal miniature model with curved blades
- 13 mm blades
- Flat handle
- Dull finish
- Overall length 100 mm

VITREORETINAL

60-549 HORIZONTAL BLADES




Squeeze Action Handle
21G / 23G

60-550 VERTICAL BLADES




Squeeze Action Handle
21G / 23G

60-551 CURVED BLADES



Squeeze Action Handle
21G / 23G

60-552 ANGLED ON SIDE




Squeeze Action Handle
21G / 23G




SCISSORS

60-822 IOL CUTTING SCISSORS



60-555 DALK SCISSOR RIGHT




60-556 DALK SCISSOR LEFT



NEEDLE HOLDERS


65-560 KALT NEEDLE HOLDER STANDARD



65-560 (C) Curved
65-560 (S) Straight

- 12 mm jaws
- Flat serrated locking handle
- Dull finish
- Overall length 140 mm


65-561 CASTROVEIJO NEEDLE HOLDER STRAIGHT



65-561 With Lock
65-562 Without Lock

- Delicate jaws
- Flat serrated handle
- Dull finish
- Overall length 140 mm


65-563 CASTROVEIJO NEEDLE HOLDER CURVED



65-563 With Lock
65-564 Without Lock

- Standard curved jaws
- Flat serrated handle
- Dull finish
- Overall length 140 mm


65-565 CASTROVEIJO NEEDLE HOLDER STRAIGHT



65-565 With Lock
65-566 Without Lock

- Standard 12mm jaws
- Flat serrated handle
- Dull finish
- Overall length 140 mm


65-567 CASTROVEIJO NEEDLE HOLDER CURVED



65-567 With Lock
65-568 Without Lock

- Standard 12mm jaws
- Flat serrated handle
- Dull finish
- Overall length 140 mm

65-569 BARRAQUER NEEDLE HOLDER STRAIGHT




65-569 With Lock
65-570 Without Lock

- Standard 12.5mm jaws
- Flat serrated handle
- Dull finish
- Overall length 125 mm

NEEDLE HOLDERS

65-571 BARRAQUER NEEDLE HOLDER, CURVED




65-571 with lock

65-572 without lock

- Standard jaws 12.5 mm
- Round knurled handle
- Overall length 125 mm

65-574 BARRAQUER NEEDLE HOLDER, STRAIGHT, 14.0 CM




65-573 with lock

65-574 without lock

- Standard jaws 12.5 mm
- Round knurled handle
- Overall length 140 mm

65-576 BARRAQUER NEEDLE HOLDER, CURVED, 14.0 CM




65-575 with lock

65-576 without lock

- Standard jaws 12.5 mm
- Round knurled handle
- Overall length 140 mm

65-580 BARRAQUER MCPHERSON NEEDLE HOLDER, CURVED




65-579 with lock

65-580 without lock

- Extra delicate jaws
- Round knurled handle
- Overall length 125 mm

65-582 BABY - BARRAQUER NEEDLE HOLDER, STRAIGHT




65-581 with lock

65-582 without lock

- Fine jaws
- Round knurled handle
- Overall length 100 mm

65-584 BABY - BARRAQUER NEEDLE HOLDER, CURVED





65-583 with lock

65-584 without lock



- Fine jaws
- Round knurled handle
- Overall length 100 mm

NEEDLE HOLDERS


65-586 PIERSE NEEDLE HOLDER, CURVED W/O LOCK




65-587 BIRKS NEEDLE HOLDER, DELICATE CURVED




65-591 JAEFFE NEEDLE HOLDER, CURVED JAWS,



65-592 JAEFFE NEEDLE HOLDER, CURVED JAWS, W/O LOCK



65-599 PIERSE NEEDLE HOLDER, CURVED





65-599 without lock

65-600 with lock

- Round knurled handle
- Overall length 140 mm

65-601 SILCOCK NEEDLE HOLDER



NEEDLE HOLDERS

BLADE BREAKER / HOLDER

70-607 SWISS MODEL RAJOR BLADE BRAKER / HOLDER 9.5CM



70-608 SWISS MODEL RAJOR BLADE BRAKER / HOLDER 12CM



70-609 BP HANDLE



STERILIZATION TRAYS

75-700 Poly Carbonate Tray With Mat / Slot

Size : 6"x2.5"x0.75"

75-702 Poly Carbonate Tray With Mat / Slot

Size : 7.5"x4"x0.75"

75-704 Poly Carbonate Tray With Mat / Slot

Size : 10"x6"x0.75"

75-707 Poly Carbonate Tray With Double Mat

Size : 10"x6"x2.25"



75-708 S.S. Tray With Mat


Size : 7.5"x5"x0.75"

75-709 S.S. Tray With Mat

Size : 10"x5"x0.75"


MARKERS

80-634 MACHATT LASIK MARKER




- Three Circle overlapping each other
- Overall length 114 mm

80-635 CHRISTENBURY LASIK MARKER




- With centering pin
- Overall length 114 mm

80-636 PERONE LASIK MARKER




- Outer Diameter 10.5 mm
- Overall length 104 mm

80-637 MALONEY LASIK FLAP FLIPPER




- 11 mm angled to tip
- 0.50 mm flat profile
- Overall length 114 mm

80-638 HOCKEY STICK KNIFE




- Used to detach epithelium from the incision by the sharp anterior part of knife and to reposit the epithelium after ablation by the blunt part of the knife

80-639 LASEK MICRO HOE




- Used to detach and lift edge of the epithelium flap

80-640 LASEK FLAP PROTECTOR




80-641 LASEK EPI BOW DISSECTOR



80-642 LASEK FLAP REPOSITOR



80-643 LASEK TREPHINE




45-643 Size 8.0 mm

45-644 Size 8.5 mm

45-645 Size 9.0 mm

• Used to create epithelium flap

80-646 LASEK ALCHOHOL WELL



45-646 Size 8.5 mm


45-647 Size 9.0 mm

45-648 Size 9.5 mm

• Used to apply alcohol solution on epithelium


MARKERS

80-650 ARNOLD MARKER




- Corneal Marker for Astigmatic Keratoplasty
- 9 mm Diameter
- Overall length 105 mm

80-651 KOCH I AXIS MARKER




- C-shaped
- Marking at 45°, 60° and 80°
- 10.5 mm Diameter
- Overall length 102 mm

80-652 KOCH II AXIS MARKER




- Corneal Marker for LRI
- Fixed marks at 90°
- 9.5 mm Diameter
- Overall length 105 mm

80-653 PRE-OP TORIC REFERENCE MARKER




80-654 FUKUYAMA LRI MARKER 60°



Marks every 10° for a total of 60° from the 9mm - 11mm diameter of the cornea

Two fixed reference point mark 90° and 270° from the steepest meridian of the cornea. Overall length 106mm, Diameter 11.5mm

80-655 FUKUYAMA LRI MARKER 120°




Marker measures 120° symmetrically on each side of the marker.

Marks every 10° for from the periphery of the cornea. Overall length 106mm, Diameter 11.5mm


MARKERS

80-661 CAPSULORHEXIS MARKER



- 5 mm Diameter
- Overall length 120 mm


80-662 KELLEN CAPSULORHEXIS MARKER



- Marking Ring With Four non-marking alignment arms.


- **80-662** 5.5 mm Diameter
- **80-663** 6.0 mm Diameter

80-664 BORES AXIS MARKER




- Overall length 120 mm

80-665 AXIS MARKER




- Overall length 100 mm

80-666 HENDERSON DEGREE GAUGE



- Degree Gauge With 2 Notches on Internal Side.
- Inside Diameter 11.7 mm
- Overall length 110 mm


80-667 FEMTO FLAP LIFTER AND RETREATMENT SPATULA



- IntraLase Flap Lifter and Retreatment Spatula


DALK-DSEK

80-668 POINTED DISSECTOR



- For lamellar Dissection
- Over all length 120 mm


80-669 TRIFACET SPATULA



- Over all length 120 mm


DSEK

80-670 CURVED DISSECTOR




- 80-670.1 angled, straight blade
- Over all length 125 mm

80-670.1 CURVED ANGLED DISSECTOR




- Over all length 125 mm

80-672 DSAEK DOUBLE ENDED MARKER 8.00/9.00 MM TITANIUM



80-674 DSAEK DEXATOME SPATULA




80-676 DSAEK DESEMENT'S STRIPPER




80-678 DSAEK STROMAL SCRUBBER



80-687 SMILE SHARP & BLUNT SPATULA




80-688 SMILE DOUBLE ENDED SPATULA





MARKERS


80-753 PRE-OP TORIC REFERENCE MARKER WITH BUBBLE LEVEL, SIDE HANDLE, TITANIUM




80-804 TWO BLADE AXIS MARKER FOR TORIC IOL, TITANIUM



80-812 AXIS MARKER WITH CENTRE PIN POINTER




80-814 GLUED IOL MARKER




80-815 INTRAVITREAL INJECTION GUIDE, TITANIUM




80-817 KAMRA 1.6MM/3.8MM DOUBLE MARKER




80-818 LENS CUTTER



80-824 DMEK BALL MARKER



80-825 DMEK "F" MARKER




80-826 DMEK "S" MARKER




80-827 DMEK "P" MARKER



80-828 STROMAL SPATULA




80-830 DMEK MICRO BALL STRIPPER



MANIPULATORS

45-271 D JAFFE-BECHERT NUCLEUS ROTATORS



- Y shaped blunt tip
- Angle 45°
- Round Disposable handle
- Overall length 130 mm

45-278 D LESTER PUSH-PULL MANIPULATOR



- 45-277 D Straight
- 45-278 D Angled

- .25 mm diameter tip
- Angle 45°
- Round Disposable handle
- Overall length 130 mm

45-285 D SINSKEY HOOK



- 45-285 D Straight
- 45-287 D Angled

- .25 mm diameter tip
- Round Disposable handle
- Overall length 130 mm

45-299 D FENZL IOL HOOK



- 45-298 D Straight
- 45-299 D Angled

- Blunt tip
- Round Disposable handle
- Overall length 130 mm

45-306 D LEWICKY LENS MANIPULATOR



- Blunt tip
- Curved Shaft
- Round Disposable handle
- Overall length 130 mm


45-314 D PHACO CHOPPER



- Blunt tip / Sharp tip
- Angle 45°
- Round Disposable handle
- Overall length 130 mm


MISCELLANEOUS INSTRUMENTS

JIR-005 MICRO LRIS RETRACTOR




- 8 mm long x .08 mm diameter
- Pack of 5

MICRO LRIS HOOK



- 1 mm wide blade
- Angle 30°
- Round Disposable Handle
- Overall length 130 mm

JES-001 EYE SHIELD



- Clear colour
- Plastic
- Pack of 10

JVS-020 SUPER SOAK EYE SPONGE



- Clear colour
- Plastic
- Pack of 5
- Overall Length 60 mm

JVS-010 CELLULOSE EYE SPONGE



- Clear colour
- Plastic
- Pack of 5
- Overall Length 60 mm

MISCELLANEOUS INSTRUMENTS

TREPHINE

JEP-001 ENDOTHELIAL DONOR CORNEAL BUTTON PUNCH



This mechanical device gives an excellent donor corneal button punched out from the endothelial side.

Features

- **Universal** - Any size of trephine from 6 mm to 11 mm can be fitted.
- **Precentred** - The trephine and the tefflon block are automatically centred for all sized trephines.
- **Spring** - After completion of the cut, the trephine retracts back Recoil System with inbuilt spring system.
- **Rust Proof** - Totally made out of non-ferrous metals. No danger of rusting even after autoclaving.

JST-0065 SUCTION TREPHINE




The outer shell of the trephine holds the recipient cornea by vacuum and the internal blade is rotated to take out the diseased button of cornea. The suction apparatus is resterilisable by dry heat, autoclaving or E.T.O. while the disposable cutting trephines are available in different sizes.

DISPOSABLE TREPHINE



**JDT-0060 to
JDT-0115**
Size : 6.00 mm to
11.00 mm

TEFLON BLOCKS



JTB-002 Size 32 mm

ARTIFICIAL IMPLANTS

ARTIFICIAL EYES



JAE-001

Sizes : Right and Left Eye
18,20,22,24,26 mm

CONFORMER



JCS-001

Sizes : Right and Left Eye
18,20,22,24,26 mm
Aperture : 1/3/5

SPHERICAL IMPLANTS



JSI-001

Sizes : 12,14,16,18,20, 22,24 MM

CASTROVEIJO IMPLANT



JCI-001

Sizes : 113,15,17,19 mm


DIAGNOSTIC STRIPS



JSS-001 Schirmer Strips
JFS-001 Fluorescein STrips
JRS-001 Rose Bengal Strips

Each Box of 100 Strips

TROCAR NEEDLE SET




23G / 25G
Set of 3 Trocar Needle with Infusion Cannula

ARTIFICIAL ANTERIOR CHAMBER



Pricon[®] Gold

The Premium Edge Blade



Pricon Gold Knives are manufactured using latest Technology to deliver consistently sharp & dimensionally true cutting edge.

Perfect Sharpness : Pricon Gold Blades are manufactured using high quality steel that is hardened and tempered to a uniform tensile strength & toughness. This process produces a structure that can be finished to ensure a consistently sharp, precise cutting edge.

Perfect Design : Pricon Gold offers unlimited geometric shapes not readily achievable with traditional manufacturing methods. These endless shape geometries are attained with unprecedented thinness and strict tolerances.

Perfect Incision : Perfect combination of sharpness and design delivers a precise and dimensionally true incision, always.

Stab Knife/Side Port

JGH - S015 Straight 15°
JGH - S022 Straight 22°
JGH - S030 Straight 30°
JGH - S045 Straight 45°




Slit Knife / Keratome (Single Bevel)

JGH - K100 1.00mm Angled 45°	JGH - K265 2.65mm Angled 45°
JGH - K120 1.20mm Angled 45°	JGH - K275 2.75mm Angled 45°
JGH - K140 1.40mm Angled 45°	JGH - K280 2.80mm Angled 45°
JGH - K160 1.60mm Angled 45°	JGH - K300 3.00mm Angled 45°
JGH - K220 2.20mm Angled 45°	JGH - K320 3.20mm Angled 45°
JGH - K240 2.40mm Angled 45°	



Clear Corneal Knife (Dual Bevel)

JGH - C160 1.60mm Angled 45°	With Depth Marking
JGH - C180 1.80mm Angled 45°	JGH - CD265 2.65mm Angled 45°
JGH - C220 2.20mm Angled 45°	JGH - CD275 2.75mm Angled 45°
JGH - C240 2.40mm Angled 45°	JGH - CD285 2.85mm Angled 45°
JGH - C265 2.65mm Angled 45°	JGH - CD300 3.00mm Angled 45°
JGH - C275 2.75mm Angled 45°	JGH - CD320 3.20mm Angled 45°
JGH - C285 2.85mm Angled 45°	
JGH - C300 3.00mm Angled 45°	
JGH - C320 3.20mm Angled 45°	



DISPOSABLE KNIVES

Implant Knife/Enlarger

JGH - B520 5.20mm Angled 45°
JGH - B550 5.50mm Angled 45°



Crescent Knife

JGH - C016 1.60mm Angled 45
JGH - C020 2.00mm Angled 45
JGH - C026 2.60mm Angled 45



MVR Knife

JGH - M019 19G Straight
JGH - M020 20G Straight
JGH - M023 23G Straight
JGH - M025 25G Straight
JGH - M019A 19G Angled
JGH - M020A 20G Angled
JGH - M023A 23G Angled
JGH - M025A 25G Angled



Paracentesis Knife

JGH - PA06 0.60mm Angled 45
JGH - PA08 0.80mm Angled 45
JGH - PA10 1.00mm Angled 45
JGH - PA11 1.10mm Angled 45



Trapezoidal Knife

JGH - 1214 1.2-1.4mm Angled 45
JGH - 1517 1.5-1.7mm Angled 45
JGH - 2730 2.7-3.0mm Angled 45



Restricted Depth Knife

JGH-0301 Restricted Depth300μ Straight
JGH-0303 Precision Depth300μ Straight
JGH-0603 Precision Depth600μ Straight



CLEAREDGE[®]

Single Use Ophthalmic Microsurgical Knife

Pricon offer microsurgical blades with advanced features like Colored handle, protective cap, light weight, sharp and consistent edges for better penetration.

Stab Knife / Side Port

For External Frown or First Scleral Incision or Making Side Port

JSH - S010 Angled 1.0 mm port
JSH - S012 Angled 1.2 mm port
JSH - S015 Straight 15°
JSH - S022 Straight 22°
JSH - S030 Straight 30°
JSH - S045 Straight 45°



Slit Knife / Keratome (Single Bevel)

Intended for making main port incision during cataract surgery

JSH - K100 1.00mm Sharp Tip Ang.45°
JSH - K120 1.20mm Sharp Tip Ang.45°
JSH - K140 1.40mm Sharp Tip Ang.45°
JSH - K160 1.60mm Sharp Tip Ang.45°
JSH - K180 1.80mm Sharp Tip Ang.45°
JSH - K220 2.20mm Sharp Tip Ang.45°
JSH - K240 2.40mm Sharp Tip Ang.45°
JSH - K265 2.65mm Sharp Tip Ang.45°
JSH - K275 2.75mm Sharp Tip Ang.45°
JSH - K220 2.20mm Sharp Tip Ang.45°
JSH - K280 2.80mm Sharp Tip Ang.45°
JSH - K300 3.00mm Sharp Tip Ang.45°
JSH - K320 3.20mm Sharp Tip Ang.45°



Clear Corneal Knife (Dual Bevel)

Intended for making clear corneal incision during cataract surgery

JSH - C120 1.20mm Sharp Tip Ang.45°
JSH - C140 1.40mm Sharp Tip Ang.45°
JSH - C160 1.60mm Sharp Tip Ang.45°
JSH - C180 1.80mm Sharp Tip Ang.45°
JSH - C220 2.20mm Sharp Tip Ang.45°
JSH - C240 2.20mm Sharp Tip Ang.45°
JSH - C265 2.65mm Sharp Tip Ang.45°
JSH - C275 2.75mm Sharp Tip Ang.45°
JSH - C285 2.85mm Sharp Tip Ang.45°
JSH - C300 3.00mm Sharp Tip Ang.45°
JSH - C320 3.20mm Sharp Tip Ang.45°



DISPOSABLE KNIVES

Implant Knife / Enlarger Knife

For I.O.L. Implant Incision or Enlarging or Slit

JSH - B350 3.5mm Blunt Tip Ang.45°

JSH - B520 5.2mm Blunt Tip Ang.45°

JSH - B550 5.5mm Blunt Tip Ang.45°



Crescent Knife

For Making Groov or Tunnel Incision

JSH - C016 Straight / Angled Bevel up

JSH - C020 Straight / Angled Bevel up

JSH - C023 Straight / Angled Bevel up

JSH - C026 Straight / Angled Bevel up



MVR Knife

For side port incision and vitreoretinal procedures

JSH - M019 19 G straight, double cutting edge

JSH - M020 20 G straight, double cutting edge

JSH - M023 23 G straight, double cutting edge

JSH - M025 25 G straight, double cutting edge

JSH - M019A 19 G angle, double cutting edge

JSH - M020A 20 G angle, double cutting edge

JSH - M023A 23 G angle, double cutting edge

JSH - M025A 25 G angle, double cutting edge



Restricted Depth Knife

Intended for making groove incision for 300 microns during cataract/keratoplasty surgery.

JSH - R350

JSH - R400

JSH - R450

JSH - R500

JSH - R600

JSH - SD300



Our manufacturing range also includes :

HYPODERMIC NEEDLES

NEEDLES BLOOD LANCET, BODY PIERCING

ALLIED SURGICAL EQUIPMENTS / OEM SUPPLIES

DISPOSABLE SYRINGES AND NEEDLES, PERFUSION SET

ALL SPECIAL TYPES OF NEEDLES LIKE : BIOPSY, L.P., EPIDURAL ETC.

Doc. No. IH/CAT/01, Rev, no.:01, Effective Date:01/12/21

ISCON SURGICALS LTD.

An ISO 9001:2015 & ISO 13485:2016 Certified Company



Works :
22/4, Heavy Industrial Area,
Jodhpur-342003 (Raj) India
Tel. : +91-291-2740553
Fax : +91-291-2740560

Head Office :
B-70, Marudhara Industrial Area,
II Phase Basni,
Jodhpur-342003 (Raj) India
Tel. : 0291-2740556-57-58
Fax.: 0291-2740560

Delhi Branch Office :
252-A, Lakshdeep Plaza, IInd Floor
Sant Nagar, East of Kailash
NEW DELHI-110065
Tel. : 0 11 32555529, 23864212
e-mail : delhi@pricon.co.in

Mumbai Branch Office :
21/202, Thakur House, K.V. Rajguru
Marg, Girgam, MUMBAI-400 004
Tel. : 0 22 32665908
e-mail : mumbai@pricon.co.in

Website : www.pricon.co.in e-mail : info@pricon.co.in